

# Aspect ISR G3 Q-CARD



**WARNING! Read operator manual prior to use**

## SETUP/PATIENT PREPARATION:

Moving table:

- 1) Unlock Casters using Caster Lock Pedal.
- 2) Position table to desired location.
- 3) Secure table mobility by locking Casters using Caster Lock Pedal.

AC Power:

- 4) Plug power cord into properly grounded outlet.
- 5) Protect power cord from damage by casters or other equipment.

Controls:

- 1) Emergency Stop Pushbutton disables all tabletop motion control when engaged.
- 2) If tabletop motion controls fail to operate, make sure that the Emergency Stop Pushbutton is NOT pushed in.

## BATTERY OPERATION (DC POWER):

- 1) The table is FULLY functional, but motions are slower.
- 2) Connect to AC power for a minimum of 8 hours daily to ensure a proper battery charge.
- 3) Failure to push in Emergency Stop Pushbutton when table is not in use and not connected to AC power for longer than 72 hours, may result in permanent damage to the batteries.



## CPR:


Return table to level, retract tabletop to minimize overhang and lower table position to a comfortable height to perform CPR.

## BRAKES:

The ISR G3 table braking system utilizes a passive magnetic brake; it remains engaged in the event of a complete loss of electrical power. Power is necessary to release brakes.

Braking force is limited by the magnet strength. Use caution when applying pressure to the side of the table as the brakes may be overcome if more than 70 pounds of force is applied.

To reduce potential motion during patient transfers with larger patients, release the table brake and move the top in the direction the patient will be pulled during loading or unloading until the table top transverse travel limit is reached. Release the panning handle to re-engage the lock and then transfer the patient.



## CLEANING AND STORAGE:

- 1) Unplug power cord and wrap it around the cord holder.
- 2) Push in Emergency Stop Pushbutton.
- 3) Set caster lock pedal to LOCK position.
- 4) Clean and disinfect table.




Scan to Access User Manuals Online

# Aspect ISR G3 Q-CARD




## TABLE TOP MOTION CONTROL PANHANDLE

**LOCK**  
PRESS DOWN ON TOP HALF OF PEDAL TO LOCK ALL MOVEMENT OF REAR CASTERS



**UNLOCK**  
PRESS DOWN ON BOTTOM HALF OF PEDAL TO ALLOW REAR CASTERS TO SWIVEL FREELY

**PANHANDLE KNOB/ TABLETOP FLOAT BRAKE**  
(Lift underside of Panhandle Knob to float table top. It will not operate if table top is not level)



**SWITCH LOCKOUT PUSHBUTTON**

(Amber light indicates that switches on pan handle are disabled)

**ELEVATE SWITCH**

(Green light to right of Tabletop Float Brake indicates table top is centered on the base laterally)

**TRENDELEBURG SWITCH**

(Green light to left of switch indicates table top is level for this motion)

**ISO-ROLL SWITCH**

(Green light to left of switch indicates table top is level for this motion)


### CASTER LOCK PEDAL POSITIONS

**EMERGENCY STOP PUSHBUTTON**



- PUSH IN to stop **ALL** tabletop motion.
- TWIST to reset.

## Table Top Motion Handset Control



It is recommended that the pendant control be disconnected during critical procedures to prevent inadvertent operation or motion, causing patient injury. The Motion Lockout Pushbutton on the Panhandle Unit **does not** disable the Handset Control.



Patients must be loaded from the side of table. Tilting or instability may occur if patient is loaded onto the pedestal end of table or the imaging end of the table.