

IDI Table Accessories Catalog

For IDI Aspect Series Tables

Quality products for surgical, interventional and diagnostic imaging




Table of Contents

Aspect ISR Table Accessories.....	3-5
Aspect 100-4T Table Accessories.....	6-8
Aspect 100-4 Table Accessories.....	9-10
Aspect 100UCPlus Urology Table Accessories.....	11-12
Aspect 100US GYN Table Accessories.....	13
Aspect 100RTL Table Accessories.....	14-15
Aspect 100RT Table Accessories.....	16-17
Aspect 100T Table Accessories.....	18-19
Aspect 100E Table Accessories.....	20-21
Aspect EZ-Access Table Accessories.....	22-23
Aspect High Tilt Table Accessories.....	23
Aspect Four-Way Float Top Tables Specifications & Accessory Matrix.....	24
Aspect Pain Management Table Specifications & Accessory Matrix.....	25
Aspect Urology & GYN Table Specifications & Accessories Matrix.....	26
IDI Customer Service Contacts.....	27

Aspect ISR G2 Table Accessories

Catalog # A100-1902



ISR table shown with X100-2492 tabletop
and optional head extension

Aspect ISR G2 Table Accessories Included:

- Tabletop pad, 2" thick
X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Rail-mounted panning hand control for 4-way float and motorized table movements. #A270-0163
- Patient Restraint Straps. # C000-0328
- Pendant Hand Control (provides secondary location for control of elevate, tilt and roll). # A100-2170

Aspect ISR G2 Table Optional Accessories:

Photo	Description	Catalog Number
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	Armboard Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007, above)	A310-056
	VAB Armboard Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244

Aspect ISR G2 Table Accessories

Photo	Description	Catalog Number
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	Headrest Extension Articulating carbon fiber headrest, Attaches to the imaging end of the ISR G2 table. Requires the table to have an X100-2492 style tabletop. Length = 12" (30.5 cm).	A100-2317
	Tabletop Catheter Tray Extension 24" wide. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-1769
	Tabletop Catheter Tray Extension 22" wide, mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-2322
	Radiation Shield Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-051
	Anesthesia Screen Holder Flexible "hoop" style. Mounts directly to A100-2317 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-0746, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2317 Headrest Extension.	C000-0965
	Anesthesia Screen Holder Rigid rod, single point mount. Adjusts vertically from 22" to 30" and rotates horizontally. (Requires one rotating side rail clamp C000-0746. If anesthesia screen is to be attached to tabletop area that has no side rail a clamp-on accessory rail is also required -	C000-0973
	Peripheral Extension Carbon fiber extension, mounts to imaging end of table only. 22" wide. Maximum weight for this extension is 25 lbs or 11.4 Kgs. - Requires tabletop X100-2492.	A100-2486

Aspect ISR G2 Table Accessories

Photo	Description	Catalog Number
	ISR Footswitch Foot control for Aspect ISRT G2 table. Foot control motions are: iso-roll, transverse travel, elevate/lower and Trendelenburg down/up.	A000-2071

Aspect ISR G2 Tabletop options

	Standard Tabletop X100-2492 19.5"/22"/24" X 80" (49.53cm/55.88cm/60.96cm x 203.20 cm). Includes short, recessed side rails 5" (13 cm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table. Catheter tray extensions are optional and sold separately, see part numbers A100-1769 and A100-2322 on page 4
	Optional Tabletop X100-2426 22"/24" X 80" (55.88cm/60.96cm x 203.20 cm). Catheter tray extension (A100-1769) may be attached to wide end of table only (pedestal end of table). Catheter tray extension is optional and sold separately. Optional headrest extension is not available with this tabletop version. Does not include recessed side rails at imaging end of table.
	Optional Tabletop available at no extra charge. Must specify with order.

Aspect 100-4T G2 Fluoroscopy Procedure Table Accessories

Catalog # A100-1416




Aspect 100-4T G2 Accessories Included:

- Tabletop pad, 2" thick
X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Patient Restraint Straps. # C000-0328
- Panning Hand Control, rail mounted. #A270-0173

Aspect 100-4T G2 Optional Accessories:

Photo	Description	Catalog Number
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	Armboard Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007, above)	A310-056
	VAB Armboard Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244


Aspect 100-4T G2 Table Accessories

Photo	Description	Catalog Number
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	Headrest Extension Articulating carbon fiber headrest, Attaches to the imaging end of the table. Requires the table to have an X100-2492 style tabletop. Length = 12" (30.5 cm).	A100-2317
	Tabletop Catheter Tray Extension 24" wide. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-1769
	Tabletop Catheter Tray Extension 22" wide, mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-2322
	Radiation Shield Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-051
	Anesthesia Screen Holder Flexible "hoop" style. Mounts directly to A100-2317 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-0746, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2317 Headrest Extension.	C000-0965
	Anesthesia Screen Holder Rigid rod, single point mount. Adjusts vertically from 22" to 30" and rotates horizontally. (Requires one rotating side rail clamp C000-0746. If anesthesia screen is to be attached to tabletop area that has no side rail a clamp-on accessory rail is also required -	C000-0973
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch circular post.	C000-0746

Aspect 100-4T G2 Table Accessories

Photo	Description	Catalog Number
	Foot control Foot control for Aspect 100-4T G2 table. Foot control motions are: elevate/lower and Trendelenburg down/up.	A100-2360
	Peripheral Extension Carbon fiber extension, mounts to imaging end of table only. 22" wide. Maximum weight for this extension is 25 lbs or 11.4 Kgs. - Requires tabletop X100-2492.	A100-2486

Aspect 100-4T G2 Tabletop Options

Photo	Description
	Standard Tabletop X100-2426 22"/24" X 80" (55.88cm/60.96cm x 203.20 cm). Catheter tray extension (A100-1769) may be attached to wide end of table only (pedestal end of table). Catheter tray extension is optional and sold separately. Optional headrest extension is not available with this tabletop version. Does not include recessed side rails at imaging end of table.
	Optional Tabletop X100-2492 19.5"/22"/24" X 80" (49.53cm/55.88cm/60.96cm x 203.20 cm). Includes short, recessed side rails 5" (13 cm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table. Catheter tray extensions are optional and sold separately, see part numbers A100-1769 and A100-2322 on page 4.
	Optional Tabletop available at no extra charge. Must specify with order.

Aspect 100-4 Fluoroscopy Procedure Table Accessories

Catalog # A100-566



Aspect 100-4 Accessories Included:

- Tabletop Pad, 2" thick. # X100-1780 for the X100-517 18" top or X100-1767 for the X100-982 23.5" top
- Patient Restraint Straps. # C000-0328
- ID410 Extender Board (Lateral arm support, with 18" wide top option only.) # A100-435
- Panning Hand Control. #A270-126

Aspect 100-4 Optional Accessories:

Photo	Description	Catalog Number
	Tabletop Catheter Tray Extension Mounts to pedestal end only.	A100-814
	Clamp-on Accessory Rail For Aspect 100-4	A100-1770
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1770, above)	A310-056
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059

Photo	Description	Catalog Number
	VAB Armboard Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244
	Radiation Shield Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-051
	ID410 Extender Board Tabletop Width Extender Board with pads: lateral arms support. Use with Aspect 100-4 tables with 18" or 22" wide tops.	18": A100-435 22": A100-780

Aspect 100-4 Tabletop options

Photo	Description
	Standard Tabletop X100-982 23.5" – 30" W x 80" L (59.5 – 76.2 x 203cm) Low attenuation, carbon fiber construction.
	Optional Tabletop X100-517 18" - 30" W x 80" L (45.7 – 76.2 x 203cm) Low attenuation, carbon fiber construction. This top comes with a removable ID410 Extender Board for adding width to the top when needed.
	Optional Tabletop available at no extra charge. Must specify with order.

Aspect 100-UCPLUS Urology & Multi-Purpose Imaging Table Catalog # A100-1517



Shown with optional leg holders
(#A100-2247).



Shown with tabletop extension
(A100-1400)

Aspect 100-UCPLUS Accessories Included:

- Tabletop Pad, 2" thick. # X100-1818
- Tabletop Extension, with pad. Extension: #A100-1400. Table extension pad: # X100-1817
- Pendant Hand Control. # K000-0158
- Foot Control. #A100-1540
- Drain Bag Support Hoop. #C000-0595
- Disposable drain bag and hose assembly, box of 20. #C000-0593
- Removable Radiation Shield Flap for perineal end of table. #Z100-1060
- Patient Restraint Straps (2). #C000-0571
- Arm board, rail mounted, with pad and strap. #A310-056
- Disposable covers for foot control, box of 50. #C000-0492
- Stirrups Mounting Rail Extender Kit (2 extensions, 2 clamps). #A100-1876



Above: Stirrups Mounting Rail Extender Kit. #A100-1876
For use with leg-holders (leg-holders are sold separately).

Aspect 100-UCPLUS Optional Accessories:

Photo	Description	Catalog Number
	Lithotomy Legholders (pair) Boot style, pneumatic assist with one-hand adjustment, with pads. For patients up to 350 lbs (158kg). Includes pair of side-rail clamps.	A100-2246
	Heavy Duty Lithotomy Legholders (pair) Heavy Duty, boot style, pneumatic assist with one-hand adjustment, with pads. For patients up to 525 lbs (238kg). Includes pair of side-rail clamps.	A100-2247
	Physician Elbow Supports (pair) Removable	C000-0594
	Armboard Rail mounted, with pad and strap 1 arm board included with table.	A310-056
	Clamp-on Accessory Rail For use with 100UC Extension. Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-2218

Aspect 100-UCPLUS Consumables:

Photo	Description	Catalog Number
	Vinyl Urology Table Covers Box of 20	C000-0645
	Uro Drain Bags with Hose Disposable drain bag and hose assembly. Box of 20	C000-0593
	Disposable Covers for Aspect 100UC Foot Controls 24" x 20" Box of 50	C000-0492
	Uro Drain Support Hoop Collapsible & removable drain bag support hoop.	C000-0595
	Uro Collection Unit 5 gallon capacity, Box of 10 10 liter capacity, Box of 10 20 liter capacity, Box of 10	5 gal: C000-0612 10 liter: C000-0905 20 liter: C000-0916

Aspect 100US GYN/Fertility Clinic Mobile Imaging Table

Catalog # A100-1479



Aspect 100US Accessories Included:

- Basic Knee Crutches #A100-2248 with pair of rotating side-rail clamps.
- Tabletop Pad, 1" thick. #X100-1687
- Foot Control. #A100-758
- Patient Restraint Straps. #C000-0328
- Detachable leg support tabletop extension with pad. Pad # X100-2351.
Weight capacity 75 lbs (34kg) #A100-2380

Aspect 100US Optional Accessories:

Photo	Description	Catalog Number
	Armboard Quick-release, rail mounted, with pad.	A310-056
	Hand control for 100US	K000-0041
	Side-rail clamp Rotating, with socket for attachment of accessories with 5/8 inch circular post.	C000-0746

Aspect 100RTL Fluoroscopic Procedure Table

Catalog # A100-1254



Aspect 100RTL Accessories Included:


- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733 or Pad #X100-2083 for Custom table top X100-1325.
- Pendant hand control, # K000-0197
- Restraint Straps #C000-0328
- Facial Cushion, crescent shape. #C000-0597
- Disposable covers for facial cushion, box of 50. #C000-0598

Aspect 100RTL Optional Accessories:

Photo	Description	Catalog Number
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059

Photo	Description	Catalog Number
	Scatter-Guard radiation shield draping For for 24" rectangular tabletop x100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096
	Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	Foot Control for 100RTL Foot Control unit controls all table motions.	A100-1961
	• Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50.	Cushion: #C000-0597 Covers: #C000-0598

Aspect 100-RTL Tabletop options

Photo	Description
	Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	Custom Tabletop X100-1325 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.

Aspect 100RT Fluoroscopic Procedure Table

Catalog # A100-870



Aspect 100RT Accessories Included:


- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733 or Pad #X100-2083 for Custom table top X100-1325.
- Pendant hand control, #K000-0195
- Patient Restraint Straps #C000-0328
- Facial Cushion, crescent shape. #C000-0597
- Disposable covers for facial cushion, box of 50. #C000-0598

Aspect 100RT Optional Accessories:

Photo	Description	Catalog Number
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	• Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50.	Cushion: #C000-0597 Covers: #C000-0598

	Description	Catalog Number
	Scatter-Guard radiation shield draping 24" rectangular tabletop x100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096
	Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	Foot Control for 100RT table Foot Control unit controls all table motions.	A100-1960

Aspect 100-RT Tabletop options

Photo	Description
	Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	Custom Tabletop X100-1325 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.

Aspect 100T Fluoroscopic Procedure Table

Catalog # A100-658




Aspect 100T Accessories Included:


- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733.
- Pendant Hand control. # K000-0054
- Patient Restraint Straps. #C000-0328
- Facial Cushion, crescent shape. # C000-0597
- Disposable covers for facial cushion, box of 50. # C000-0598

Aspect 100T Optional Accessories:

Photo	Description	Catalog Number
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	Scatter-Guard radiation shield draping For 24" rectangular tabletop x100-995. Dimensions: 59" long. 24" Draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096

Photo	Description	Catalog Number
	Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	Foot Control for 100T table (Substitute for Hand Control, cannot use hand and foot control at the same time.)	A100-500
	• Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50.	Cushion: #C000-0597 Covers: #C000-0598

Aspect 100T Tabletop options

Photo	Description
	Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	Custom Tabletop X100-1325 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.

Aspect 100E Fluoroscopic Procedure Table

Catalog # A100-402



Aspect 100E Accessories Included:

- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733.
- Pendant Hand control. # K000-0041
- Patient Restraint Straps. #C000-0328
- Facial Cushion, crescent shape. # C000-0597
- Disposable covers for facial cushion, box of 50. # C000-0598

Aspect 100E Optional Accessories:

Photo	Description	Catalog Number
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	Scatter-Guard radiation shield draping For 24" rectangular tabletop x100-995. Dimensions: 59" long.24" Draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096

Photo	Description	Catalog Number
	Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	Foot Control (Substitute for Hand Control, cannot use pendant hand control and foot control at the same time.)	A100-758
	• Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50.	Cushion: #C000-0597 Covers: #C000-0598

Aspect 100E Tabletop options

Photo	Description
	Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	Custom Tabletop X100-1325 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.

Aspect EZ-Access Fluoroscopic Procedure Table

Catalog # A100-2100



Shown with Facial Cutout tabletop X100-733

Aspect EZ-Access Accessories Included:


- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733.
- Pendant Hand control. # K000-0041
- Patient Restraint Straps. #C000-0328
- Facial Cushion, crescent shape. # C000-0597
- Disposable covers for facial cushion, box of 50. # C000-0598

Aspect EZ-Access Optional Accessories:

Photo	Description	Catalog Number
	Clamp-on Accessory Rail Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	Armboard Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	Foot Control (Substitute for Hand Control, cannot use pendant hand control and foot control at the same time.)	A100-758
	Scatter-Guard radiation shield draping For 24" rectangular tabletop x100-995. Dimensions: 59" long. 24" Draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096

	Scatter-Guard radiation shield draping for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table.	A610-092
---	--	----------

Aspect 100EZ Tabletop options

Photo	Description
	Standard Tabletop X100-995 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	Optional Tabletop X100-733 with facial cutout 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	Custom Tabletop X100-1325 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.

Aspect High-Tilt™ Fluoroscopic Procedure Table

Motorized 75° Upright Tilt, Fixed Height, with Built-in Footboard for Step-On Patient Loading

Catalog # A100-1708



Aspect High-Tilt™ Accessories Included:

- Tabletop pad, 2" thick with hood. #X100-2008
- Patient handgrips. #A100-171
- Restraint straps. #C000-0571
- Foot switch control. #A100-1823

Aspect Four-Way Float Table Specifications & Accessories

•=Standard feature ◇=Optional NA=Not available	100 - 4	100 - 4T G2	ISR G2
Elevate (Height without pad)	27.9 - 37.7" 71 - 96 cm	29 - 38.9" 73.7 - 99 cm	30 - 38" 76 - 96.5 cm
Trendelenburg Tilt	NA	±12°	±12°
Lateral Tilt	NA	NA	±12°
Longitudinal Tabletop Motion	40" float 102 cm	32" float 81 cm	32" float 81 cm
Transverse Tabletop Motion	10" float 25 cm	8" float 20 cm	8" float 20 cm
Standard Tabletop Size	18/30" x 80" 46/76 x 203 cm	24" x 84" 61 x 213 cm	22"/24" x 84" 56/61 x 213 cm
Optional Tabletop Size	23.5-30" x 80" 59-76 x 203 cm	NA	NA
Imaging Area	72" 182.9 cm	62" 157.5 cm	61" (154.9 cm) x 73" (185.4 cm) with headrest attached
Patient Capacity	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg
Table Weight	490 lb. 222 kg	605 lb. 274 kg	605 lb. 274 kg
Tabletop Pad (2" thick)	•	•	•
Restraint Straps	•	•	•
Pendant Hand Control	NA	NA	•
4-way Panning Hand Control	•	•	•
Foot Control	NA	(Custom order)	(Custom order)
Armboard, rail mounted	NA	◇	◇
Armboard, carbon fiber	◇	◇	◇
VAB Armboard	◇ Use of the VAB arm board with the 100-4 is only recommended on the models with 23.5" the table top . Use is only possible with a patient's weight securing the VAB to the table.		
Clamp-on Accessory Rail	◇	◇	◇
Catheter Tray extension	◇	◇	◇
Headrest extension	NA	◇	◇
Table-Mounted Radiation Shield	◇	◇	◇
Anesthesia Screen Holder	NA	◇	◇

Aspect Pain Management Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	EZ-Access	100E	100T	100RT	100RTL
Elevate (Height without pad)	19 – 38" 48 – 96.5 cm	26 - 37.8" 66 - 96 cm	28 - 39.8" 71 - 101 cm	30.4 - 42.2" 77 - 107 cm	31.2 - 41.2" 79 – 104.6 cm
Trendelenburg Tilt	NA	NA	±15°	±15°	±15°
Lateral Tilt	NA	NA	NA	20°	20°
Longitudinal Tabletop Motion	NA	NA	NA	NA	10"
Standard Tabletop Size	24" x 80" 61 x 203cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 – 203 cm
Optional Tabletop Size	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 – 203 cm
Optional Tabletop Size	NA	NA	NA	28" x 80' 72 x 203 cm	28" X 80" 72 x 203 cm
Imaging Area	68" 172.7 cm	66" 167.6 cm	62" 157.5 cm	62" 157.5 cm	59" 149.9 cm
Patient Capacity	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lbs. 204 kg
Table Weight	490 lb. 222 kg	290 lb. 131 kg	325 lb. 147 kg	425 lb. 192 kg	355 lbs. 161 kg
Tabletop Pad (2" thick)	•	•	•	•	•
Restraint Straps	•	•	•	•	•
Facial Support Cushion	•	•	•	•	•
Pendant Hand Control	•	•	•	•	•
Foot Control	◊	◊	◊	◊	◊
Armboard, rail mounted	◊	◊	◊	◊	◊
Armboard, carbon fiber	◊	◊	◊	◊	◊
VAB Armboard	◊	◊	◊	◊	◊
Clamp-on Accessory Rail	◊	◊	◊	◊	◊
Scatter-Guard Radiation Shield	◊	◊	◊	◊	◊

Aspect Urology & GYN Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	100 UC-Plus	100 US
Elevate (Height without pad)	31.3" – 41.1" (79.4 – 104.4 cm)	26" – 38" (66 – 96.5cm)
Trendelenburg Tilt	± 15°	NA
Transverse travel	± 3.5", 7" total (± 8.9 cm, 17.8 cm total)	NA
Longitudinal Tabletop Motion	10" (25.4 cm)	NA
Tabletop Size	28"W x 46"L (71 x 117 cm)	24"W x 46"L (61 x 117 cm)
Tabletop Extension Size	24"W x 34"L (61 x 86 cm)	NA
Imaging Area	Main tabletop: 25"W x 20"L (63.5 x 50.8 cm)	22"W x 32"L (56 x 81 cm)
Imaging Area with table extension	54" (137 cm)	NA
Patient Weight Capacity	Upper body on main table: 525lb. (238 kg) Upper body on extension: 400 lbs (181kg) Load capacity on extension only: 200 lbs (90.7kg)	450 lbs (204 kg)
Table Weight	490 lb. (222 kg)	375 lbs ((170kg)
Tabletop Pad (2" thick)	•	•
Tabletop extension with pad	•	Leg support tabletop extension
Hand control pendant & footswitch control	•	Footswitch only
Patient restraint straps	•	•
Removeable drain bag support hoop	•	◊
Disposable drain bag & hose assembly, box of 20	•	◊
Removable radiation shield flap for perineal end of table	•	◊
Armboard, rail-mounted with pad and strap	•	◊
Disposable covers for footswitch, box of 50	•	◊
Radiographic cassette tray	•	NA

IDI Product Support & Customer Service
Hours: 8:00 AM-5:00 PM Eastern Standard Time
Phone: 978-829-0009
Toll Free: 877-304-5434
Fax: 978-829-0027
Email: customerservice@imagediagnostics.com
ImageDiagnostics.com

Image Diagnostics, Inc.

Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

Image Diagnostics, Inc.

310 Authority Drive
Fitchburg, Massachusetts
01420

Toll free: 877.304.5434
Main: 978.829.0009
Fax: 978.829.0027

ImageDiagnostics.com

Please direct all inquiries to:
sales@imagediagnostics.com

idi
Image Diagnostics Inc.



A Veteran-Owned Business