

# IDI Table Accessories Catalog

For IDI Aspect Series Tables

Quality products for surgical, interventional and diagnostic imaging



## Table of Contents

Product/Service	Page
Aspect ISR G2 Table Accessories	3-5
Aspect 100-4T G2 Table Accessories	6-8
Aspect 100-4 Table Accessories	9-10
Aspect 100UC <i>Plus</i> Urology Table Accessories	11-13
Aspect 100US GYN Table Accessories	14
Aspect 100RTL Table Accessories	15-16
Aspect 100RT Table Accessories	17-18
Aspect 100T Table Accessories	19-20
Aspect 100E Table Accessories	21-22
Aspect High Tilt Table Accessories	23
Aspect Four-Way Float Top Tables Specifications & Accessory Matrix	24
Aspect Pain Management Table Specifications & Accessory Matrix	25
Aspect Urology & GYN Table Specifications & Accessory Matrix	26
IDI Customer Service Contacts	27

## Aspect ISR G2 Table Accessories

### Catalog # A100-1902






ISR table shown with X100-2492 tabletop  
and optional head extension

### Aspect ISR G2 Table Accessories Included:

- Tabletop pad, 2" thick  
# X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Rail-mounted panning hand control for 4-way float and motorized table movements. #A270-0163
- Patient Restraint Straps. # C000-0328
- Pendant Hand Control (provides secondary control of elevate, trendelenburg and roll). # A100-2170

### Aspect ISR G2 Table Optional Accessories:


Photo	Description	Catalog Number
	<b>Clamp-on Accessory Rail</b> Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	<b>Armboard</b> Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007, above)	A310-056
	<b>VAB Armboard</b> Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244

## Aspect ISR G2 Table Accessories



Photo	Description	Catalog Number
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	<b>Headrest Extension</b> Articulating carbon fiber headrest, Attaches to the imaging end of the ISR G2 table. Requires the table to have an X100-2492 style tabletop. Length = 12" (30.5 cm).	A100-2317
	<b>Tabletop Catheter Tray Extension</b> 24" wide. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-1769
	<b>Tabletop Catheter Tray Extension</b> 22" wide, mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-2322
	<b>Radiation Shield</b> Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-051
	<b>Anesthesia Screen Holder</b> Flexible "hoop" style. Mounts directly to A100-2317 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-0746, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2317 Headrest Extension.	C000-0965
	<b>Anesthesia Screen Holder</b> Rigid rod, single point mount. Adjusts vertically from 22" to 30" and rotates horizontally. (Requires one rotating side rail clamp C000-0746. If anesthesia screen is to be attached to tabletop area that has no side rail a clamp-on accessory rail is also required -	C000-0973
	<b>Peripheral Extension</b> Carbon fiber extension, mounts to imaging end of table only. 22" wide. Maximum weight for this extension is 25 lbs or 11.4 Kgs. - Requires tabletop X100-2492.	A100-2486



## Aspect ISR G2 Table Accessories

Photo	Description	Catalog Number
	<b>ISR Table Foot Switch with Safety Hood.</b> Foot control for Aspect ISR G2 table. Foot control motions are: iso-roll, transverse travel, elevate/lower and Trendelenburg down/up. (Safety hood helps avoid accidental motorized table movement from foot switch).	A000-2071

## Aspect ISR G2 Tabletop options

	<b>Standard Tabletop X100-2492</b> 19.5"/22"/24" X 84" (49.5cm/55.8cm/60.96cm x 213.3 cm). Includes short, recessed side rails 5" (13 cm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table. Catheter tray extensions are optional and sold separately, see part numbers A100-1769 and A100-2322 on page 4
	<b>Optional Tabletop X100-2426</b> 22"/24" X 84" (55.8cm/60.9cm x 213.3 cm). Catheter tray extension (A100-1769) may be attached to wide end of table only (pedestal end of table). Catheter tray extension is optional and sold separately. Optional headrest extension is not available with this tabletop version. <b>Does not</b> include recessed side rails at imaging end of table.
	Optional Tabletop available at no extra charge. Must specify with order.

## Aspect 100-4T G2 Fluoroscopy Procedure Table Accessories




### Catalog # A100-1416



### Aspect 100-4T G2 Accessories Included:

- Tabletop pad, 2" thick  
# X100-2310 for the X100-2492 table top or # X100-2489 for the X100-2426 table top.
- Patient Restraint Straps. # C000-0328
- Panning Hand Control, rail mounted. #A270-0173



### Aspect 100-4T G2 Optional Accessories:

Photo	Description	Catalog Number
	<b>Clamp-on Accessory Rail</b> Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	<b>Armboard</b> Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007, above)	A310-056
	<b>VAB Armboard</b> Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244



## Aspect 100-4T G2 Table Accessories

Photo	Description	Catalog Number
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	<b>Headrest Extension</b> Articulating carbon fiber headrest, Attaches to the imaging end of the table. Requires the table to have an X100-2492 style tabletop. Length = 12" (30.5 cm).	A100-2317
	<b>Tabletop Catheter Tray Extension</b> 24" wide. Mounts to pedestal end of table only. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-1769
	<b>Tabletop Catheter Tray Extension</b> 22" wide, mounts to imaging end of table only. Requires tabletop X100-2492. The maximum weight this extension can support is 20 LBS or 9 Kgs.	A100-2322
	<b>Radiation Shield</b> Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-051
	<b>Anesthesia Screen Holder</b> Flexible "hoop" style. Mounts directly to A100-2317 Headrest Extension. For use on accessory rails - a pair of rotating side rail clamps, part number C000-0746, are required. Rotating side rail clamps and clamp-on accessory rails are not required if table is ordered with the A100-2317 Headrest Extension.	C000-0965
	<b>Anesthesia Screen Holder</b> Rigid rod, single point mount. Adjusts vertically from 22" to 30" and rotates horizontally. (Requires one rotating side rail clamp C000-0746. If anesthesia screen is to be attached to tabletop area that has no side rail a clamp-on accessory rail is also required -	C000-0973
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 5/8 inch circular post.	C000-0746

## Aspect 100-4T G2 Table Accessories

Photo	Description	Catalog Number
	<b>Foot control</b> Foot control for Aspect 100-4T G2 table. Foot control motions are: elevate/lower and Trendelenburg down/up.	A100-2360
	<b>Peripheral Extension</b> Carbon fiber extension, mounts to imaging end of table only. 22" wide. Maximum weight for this extension is 25 lbs or 11.4 Kgs. - Requires tabletop X100-2492.	A100-2486

## Aspect 100-4T G2 Tabletop Options

Photo	Description
	<b>Standard Tabletop X100-2492</b> 19.5"/22"/24" X 84" (49.5cm/55.8cm/60.96cm x 213.3 cm). Includes short, recessed side rails 5" (13 cm) in length at imaging end of table that facilitate use of removable headrest extension (option) and also facilitates attachment of catheter tray extension at narrow end of table. Two separate catheter tray extensions are available: One for the pedestal end of table and one for the imaging end of table. Catheter tray extensions are optional and sold separately, see part numbers A100-1769 and A100-2322 on page 4.
	<b>Optional Tabletop X100-2426</b> 22"/24" X 84" (55.8cm/60.9cm x 213.3 cm). Catheter tray extension (A100-1769) may be attached to wide end of table only (pedestal end of table). Catheter tray extension is optional and sold separately. Optional headrest extension is not available with this tabletop version. <b>Does not</b> include recessed side rails at imaging end of table.
	Optional Tabletop available at no extra charge. Must specify with order.

## Aspect 100-4 Fluoroscopy Procedure Table Accessories

### Catalog # A100-566



#### Aspect 100-4 Accessories Included:

- Tabletop Pad, 2" thick. # X100-1780 for the X100-517 18" top or X100-1767 for the X100-982 23.5" top
- Patient Restraint Straps. # C000-0328
- ID410 Extender Board (Lateral arm support, with 18" wide top option only.) # A100-435
- Panning Hand Control. #A270-126

#### Aspect 100-4 Optional Accessories:









Photo	Description	Catalog Number
	<b>Tabletop Catheter Tray Extension</b> Mounts to pedestal end only.	A100-814
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007, above)	A310-056
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	<b>VAB Armboard</b> Carbon fiber, one-piece, extra wide armboard for vascular access procedures (body section slides under patient/table pad). Includes strap and pad.	A100-2244

Photo	Description	Catalog Number
	<b>Radiation Shield</b> Articulating, rail-mounted. Attaches to table side-rail at base end of table, either side, adjustable. Articulating joint allows shield to conform to table contour. Snap-on upper section above level of tabletop is removable.	A610-051
	<b>ID410 Extender Board</b> Tabletop Width Extender Board with pads: lateral arms support. Use with Aspect 100-4 tables with 18" or 22" wide tops.	18": A100-435 22": A100-780

### Aspect 100-4 Tabletop options


Photo	Description
	<b>Standard Tabletop X100-982</b> 23.5" – 30" W x 80" L (59.5 – 76.2 x 203cm) Low attenuation, carbon fiber construction.
	<b>Optional Tabletop X100-517</b> 18" - 30" W x 80" L (45.7 – 76.2 x 203cm) Low attenuation, carbon fiber construction. This top comes with a removable ID410 Extender Board for adding width to the top when needed.
	Optional Tabletop available at no extra charge. Must specify with order.



## Aspect 100-UCPLUS Urology & Multi-Purpose Imaging Table Catalog # A100-1517





Shown with optional leg holders  
(#A100-2247).



Shown with tabletop extension  
(A100-1400)

### Aspect 100-UCPLUS Accessories Included:

- Tabletop Pad, 2" thick. # X100-1818
- Tabletop Extension, with pad. Extension: #A100-1400. Table extension pad: # X100-1817
- Pendant Hand Control. # K000-0158
- Foot Control. #A100-1540
- Drain Bag Support Hoop. #C000-0595
- Disposable drain bag and hose assembly, box of 20. #C000-0593
- Removable Radiation Shield Flap for perineal end of table. #Z100-1060
- Patient Restraint Straps (2). #C000-0571
- Arm board, rail mounted, with pad and strap. #A310-056
- Disposable covers for foot control, box of 50. #C000-0492
- Stirrups Mounting Rail Extender Kit (2 extensions, 2 clamps). #A100-1876








Above: Stirrups Mounting Rail Extender Kit. #A100-1876  
For use with leg-holders (leg-holders are sold separately).

## Aspect 100-UCPLUS Optional Accessories:

Photo	Description	Catalog Number
	<b>Lithotomy Legholders (pair)</b> Boot style, pneumatic assist with one-hand adjustment, with pads. For patients up to 350 lbs (158kg). Includes pair of side-rail clamps.	A100-2246
	<b>Heavy Duty Lithotomy Legholders (pair)</b> Heavy Duty, boot style, pneumatic assist with one-hand adjustment, with pads. For patients up to 525 lbs (238kg). Includes pair of side-rail clamps.	A100-2247
	<b>Physician Elbow Supports (pair)</b> Removable	C000-0594
	<b>Armboard</b> Rail mounted, with pad and strap 1 arm board included with table.	A310-056
	<b>Clamp-on Accessory Rail</b> For use with 100UC Extension. Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-2218
	<b>Split Leg Holders</b> For use with Aspect 100UC Plus Tables. Includes pads. Prone percutaneous, lateral or modified lateral nephrostomy. Includes two split leg holder clamps # C000-1115.	C000-1146
	<b>Replacement Split Leg Holder Clamp</b> For use with split leg holders. # C000-1146	C000-1115

## Aspect 100-UCPLUS Consumables:

Photo	Description	Catalog Number
	<b>Vinyl Urology Table Covers</b> Box of 20	C000-0645
	<b>Uro Drain Bags with Hose</b> Disposable drain bag and hose assembly. Box of 20	C000-0593
	<b>Disposable Covers for Aspect 100UC Foot Controls</b> 24" x 20" Box of 50	C000-0492
	<b>Uro Drain Support Hoop</b> Collapsible & removable drain bag support hoop.	C000-0595
	<b>Uro Collection Unit</b> 5 gallon capacity, Box of 10 10 liter capacity, Box of 10 20 liter capacity, Box of 10	5 gal: C000-0612 10 liter: C000-0905 20 liter: C000-0916

## Aspect 100US GYN/Fertility Clinic Mobile Imaging Table




Catalog # A100-1479



### Aspect 100US Accessories Included:

- Basic Knee Crutches #A100-2248 with pair of rotating side-rail clamps.
  - Tabletop Pad, 1" thick. #X100-1687
  - Foot Control. #A100-758
  - Patient Restraint Straps. #C000-0328
  - Detachable leg support tabletop extension with pad. # X100-2351.
- Weight capacity 75 lbs (34kg)

### Aspect 100US Optional Accessories:

Photo	Description	Catalog Number
	<b>Armboard</b> Quick-release, rail mounted, with pad.	A310-056
	<b>Hand control for 100US</b>	K000-0041
	<b>Side-rail clamp</b> Rotating, with socket for attachment of accessories with 5/8 inch circular post.	C000-0746

## Aspect 100RTL Fluoroscopic Procedure Table

### Catalog # A100-1254



### Aspect 100RTL Accessories Included:

- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733 or Pad #X100-2083 for Custom table top X100-1325.
- Pendant hand control, # K000-0197
- Restraint Straps #C000-0328
- Facial Cushion, crescent shape. #C000-0597
- Disposable covers for facial cushion, box of 50. #C000-0598

### Aspect 100RTL Optional Accessories:











Photo	Description	Catalog Number
	<b>Clamp-on Accessory Rail</b> Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059

Photo	Description	Catalog Number
	<b>Scatter-Guard radiation shield draping</b> For for 24" rectangular tabletop x100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096
	<b>Scatter-Guard radiation shield draping</b> for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	<b>Foot Control for 100RTL</b> Foot Control unit controls all table motions.	A100-1961
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	Cushion: #C000-0597  Covers: #C000-0598

### Aspect 100-RTL Tabletop options

Photo	Description
	<b>Standard Tabletop X100-995</b> 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	<b>Optional Tabletop X100-733 with facial cutout</b> 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	<b>Custom Tabletop X100-1325</b> 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.



## Aspect 100RT Fluoroscopic Procedure Table





### Catalog # A100-870



### Aspect 100RT Accessories Included:




- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733 or Pad #X100-2083 for Custom table top X100-1325.
- Pendant hand control, #K000-0195
- Patient Restraint Straps #C000-0328
- Facial Cushion, crescent shape. #C000-0597
- Disposable covers for facial cushion, box of 50. #C000-0598

### Aspect 100RT Optional Accessories:

Photo	Description	Catalog Number
	<b>Clamp-on Accessory Rail</b> Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	Cushion: #C000-0597  Covers: #C000-0598

	Description	Catalog Number
	<b>Scatter-Guard radiation shield draping</b> 24" rectangular tabletop x100-995. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096
	<b>Scatter-Guard radiation shield draping</b> for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	<b>Foot Control for 100RT table</b> Foot Control unit controls all table motions.	A100-1960

## Aspect 100-RT Tabletop options

Photo	Description
	<b>Standard Tabletop X100-995</b> 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	<b>Optional Tabletop X100-733 with facial cutout</b> 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	<b>Custom Tabletop X100-1325</b> 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.

## Aspect 100T Fluoroscopic Procedure Table

### Catalog # A100-658



### Aspect 100T Accessories Included:

- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733.
- Pendant Hand control. # K000-0054
- Patient Restraint Straps. #C000-0328
- Facial Cushion, crescent shape. # C000-0597
- Disposable covers for facial cushion, box of 50. # C000-0598

### Aspect 100T Optional Accessories:











Photo	Description	Catalog Number
	<b>Clamp-on Accessory Rail</b> Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	<b>Scatter-Guard radiation shield draping</b> For 24" rectangular tabletop x100-995. Dimensions: 59" long. 24" Draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096

Photo	Description	Catalog Number
	<b>Scatter-Guard radiation shield draping</b> for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	<b>Foot Control for 100T table</b> (Substitute for Hand Control, cannot use hand and foot control at the same time.)	A100-500
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	Cushion: #C000-0597  Covers: #C000-0598

## Aspect 100T Tabletop options

Photo	Description
	<b>Standard Tabletop X100-995</b> 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	<b>Optional Tabletop X100-733 with facial cutout</b> 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	<b>Custom Tabletop X100-1325</b> 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.

## Aspect 100E Fluoroscopic Procedure Table

### Catalog # A100-402






### Aspect 100E Accessories Included:




- Tabletop Pad, 2" thick. Pad #X100-1742 for standard table top X100-995 or Pad #X100-1757 for optional table top X100-733.
- Pendant Hand control. # K000-0041
- Patient Restraint Straps. #C000-0328
- Facial Cushion, crescent shape. # C000-0597
- Disposable covers for facial cushion, box of 50. # C000-0598

### Aspect 100E Optional Accessories:

Photo	Description	Catalog Number
	<b>Clamp-on Accessory Rail</b> Rail is 6" long x 3/8" thick x 1-1/8" high.	A100-1007
	<b>Armboard</b> Quick-release, rail mounted. Requires Clamp-on Accessory Rail A100-1007.	A310-056
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient/table pad), flat. Includes pad and strap.	A310-059
	<b>Scatter-Guard radiation shield draping</b> For 24" rectangular tabletop x100-995. Dimensions: 59" long.24" Draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop <i>Note: Not for use with arm board A310-056.</i>	A610-096

Photo	Description	Catalog Number
	<b>Scatter-Guard radiation shield draping</b> for 22" tabletop X100-733 models with 22" tapered-end tabletop with facial cutout. Dimensions: 59" long. 24" draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop. <i>Note: Not for use with arm board A310-056.</i>	A610-092
	<b>Foot Control</b> (Substitute for Hand Control, cannot use pendant hand control and foot control at the same time.)	A100-758
	<ul style="list-style-type: none"> <li>• Facial Cushion, crescent shape.</li> <li>• Disposable covers for facial cushion, box of 50.</li> </ul>	Cushion: #C000-0597  Covers: #C000-0598

## Aspect 100E Tabletop options

Photo	Description
	<b>Standard Tabletop X100-995</b> 24" x 80" (61 x 203cm) Low attenuation, carbon fiber construction. Standard tops may be substituted with optional sizes/styles as shown below at no extra charge.
	<b>Optional Tabletop X100-733 with facial cutout</b> 22" x 80" (56 x 203cm). Low attenuation, carbon fiber construction. Optional tabletop available at no extra charge. Must specify top selection with order.
	<b>Custom Tabletop X100-1325</b> 28" x 80" (71 x 203cm) Low attenuation, carbon fiber construction. This custom top is a special order item. Pricing on request.



## Aspect High-Tilt™ Fluoroscopic Procedure Table

Motorized 75° Upright Tilt, Fixed Height, with Built-in Footboard for Step-On Patient Loading

Catalog # A100-1708



### Aspect High-Tilt™ Accessories Included:

- Tabletop pad, 2" thick with hood. #X100-2008
- Patient handgrips. #A100-1711
- Restraint straps. #C000-0571
- Foot switch control. #A100-1823

## Aspect Four-Way Float Table Specifications & Accessories

•=Standard feature ◊=Optional NA=Not available	100 - 4	100 - 4T G2	ISR G2
<b>Elevate (Height with 2" pad)</b>	29.9 - 39.7" 75.9 - 100.8 cm	32 - 42" 81 - 107 cm	32 - 40" 81 - 102 cm
<b>Trendelenburg Tilt</b>	NA	±12°	±12°
<b>Lateral Tilt</b>	NA	NA	±12°
<b>Longitudinal Tabletop Motion</b>	40" float 102 cm	32" float 81 cm	32" float 81 cm
<b>Transverse Tabletop Motion</b>	10" float 25 cm	8" float 20 cm	8" float 20 cm
<b>Standard Tabletop Size</b>	18/30" x 80" 46/76 x 203 cm	24"/22"/19.5" x 84" 60/55.8/49.5 x 213 cm	24"/22"/19.5" x 84" 60/55.8/49.5 x 213 cm
<b>Optional Tabletop Size</b>	23.5-30" x 80" 59-76 x 203 cm	22"/24" x 84" 56 /61 x 213 cm	22"/24" x 84" 56 /61 x 213 cm
<b>Imaging Area</b>	72" 182.9 cm	61" 155 cm	61" (155cm), 74" (185.4 cm) with optional peripheral extension
<b>Patient Capacity</b>	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg
<b>Table Weight</b>	490 lb. 222 kg	605 lb. 274 kg	605 lb. 274 kg
<b>Tabletop Pad (2" thick)</b>	•	•	•
<b>Restraint Straps</b>	•	•	•
<b>Pendant Hand Control</b>	NA	NA	•
<b>4-way Panning Hand Control</b>	•	•	•
<b>Foot Control</b>	NA	(Custom order)	(Custom order)
<b>Armboard, rail mounted</b>	NA	◊	◊
<b>Armboard, carbon fiber</b>	◊	◊	◊
<b>VAB Armboard</b>	◊ Use of the VAB arm board with the 100-4 is only recommended on the models with 23.5" the table top . Use is only possible with a patient's weight securing the VAB to the table.		
<b>Clamp-on Accessory Rail</b>	◊	◊	◊
<b>Catheter Tray extension</b>	◊	◊	◊
<b>Headrest extension</b>	NA	◊	◊
<b>Table-Mounted Radiation Shield</b>	◊	◊	◊
<b>Anesthesia Screen Holder</b>	NA	◊	◊

## Aspect Pain Management Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	100E	100T	100RT	100RTL
Elevate (Height with 2" pad)	28 - 39.8" 71 - 101 cm	30 - 41.8" 76 - 106 cm	32.4 - 44.2" 82 - 112 cm	31.2 - 41.2" 79.2 - 104.6 cm
Trendelenburg Tilt	NA	±15°	±15°	±15°
Lateral Tilt	NA	NA	20°	20°
Longitudinal Tabletop Motion	NA	NA	NA	10"
Standard Tabletop Size	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 - 203 cm
Optional Tabletop Size	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 - 203 cm
Custom Tabletop Size	28" x 80" 72 x 203 cm	28" x 80" 72 x 203 cm	28" x 80" 72 x 203 cm	28" X 80" 72 x 203 cm
Imaging Area	66" 167.6 cm	62" 157.5 cm	62" 157.5 cm	59" 149.9 cm
Patient Capacity	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lbs. 204 kg
Table Weight	290 lb. 131 kg	325 lb. 147 kg	425 lb. 192 kg	355 lbs. 161 kg
Tabletop Pad (2" thick)	•	•	•	•
Restraint Straps	•	•	•	•
Facial Support Cushion	•	•	•	•
Pendant Hand Control	•	•	•	•
Foot Control	◊	◊	◊	◊
Armboard, rail mounted	◊	◊	◊	◊
Armboard, carbon fiber	◊	◊	◊	◊
VAB Armboard	◊	◊	◊	◊
Clamp-on Accessory Rail	◊	◊	◊	◊
Scatter-Guard Radiation Shield	◊	◊	◊	◊

## Aspect Urology & GYN Table Specifications & Accessories

•Standard feature ◊Optional NA Not available	100 UC-Plus	100 US
Elevate (Height with pad)	33.3" – 43.1" (84.5 – 109.4 cm)	27" – 39" (68.5 – 99cm)
Trendelenburg Tilt	± 15°	NA
Transverse travel	± 3.5", 7" total (± 8.9 cm, 17.8 cm total)	NA
Longitudinal Tabletop Motion	10" (25.4 cm)	NA
Tabletop Size	28"W x 46"L (71 x 117 cm)	24"W x 46"L (61 x 117 cm)
Tabletop Extension Size	24"W x 34"L (61 x 86 cm)	NA
Imaging Area	Main tabletop: 25"W x 20"L ( 63.5 x 50.8 cm)	22"W x 32"L (56 x 81 cm)
Imaging Area with table extension	54" (137 cm)	NA
Patient Weight Capacity	Upper body on main table: 525lb. (238 kg) Upper body on extension: 400 lbs ( 181kg) Load capacity on extension only: 200 lbs (90.7kg)	450 lbs (204 kg)
Table Weight	490 lb. (222 kg)	375 lbs ((170kg)
Tabletop Pad (2" thick)	•	•
Tabletop extension with pad	•	Leg support tabletop extension
Hand control pendant & footswitch control	•	Footswitch only
Patient restraint straps	•	•
Removeable drain bag support hoop	•	◊
Disposable drain bag & hose assembly, box of 20	•	◊
Removable radiation shield flap for perineal end of table	•	◊
Armboard, rail-mounted with pad and strap	•	◊
Split Leg Stirrups	•	NA
Disposable covers for footswitch, box of 50	•	◊
Radiographic cassette tray	•	NA

## Image **Diagnostics, Inc.**

Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

### **Image Diagnostics, Inc.**

310 Authority Drive  
Fitchburg, Massachusetts  
01420

Toll free: 877.304.5434

Main: 978.829.0009

Fax: 978.829.0027

[ImageDiagnostics.com](http://ImageDiagnostics.com)

Please direct all inquiries to:  
[sales@imagediagnostics.com](mailto:sales@imagediagnostics.com)

# idi

Image **Diagnostics Inc.**



A Veteran-Owned Business