



Product Data Sheet  
**Aspect 100US Gyn/Fertility Clinic Mobile Imaging Table**

*Ideal C-arm imaging solution for hysterosalpingograms*



- Carbon fiber tabletop with low x-ray attenuation for superior fluoroscopic imaging
- Unobstructed fluoroscopy imaging area: 22"W x 32"L (56 x 81 cm)
- Motorized tabletop elevation with footswitch control: 26" – 38" (66 – 96.5cm)
- 24"W x 46"L tabletop (61 x 117 cm)
- Patient weight capacity: 450 lbs. (204 kg)
- Standard US table side-rail: compatible with many stirrups and stirrup clamps.
- 120 VAC operation with battery back-up operation (230VAC optional)
- Easy table transport with 5" locking casters. Table weight = 375 lbs. (170 kg)
- 3-year parts warranty
- Accessories included:
  - tabletop pad
  - patient restraint straps
  - footswitch control
- Optional accessories available:
  - Lithotomy legholders (boot style or under-calf style) and attaching clamps
  - Table hand control
  - Armboards

D000-517A