


VAB

Radiolucent Arm board

Catalog # A100-2244

For vascular access and angiographic imaging of the arm.
For use only with IDI imaging tables with either 22" or 24" wide table tops.



Ultra-rigid carbon fiber construction

- High-strength, one-piece design.
- Metal-free and low x-ray attenuation.
- Arm board is held securely in place by attachment to the opposite side of the tabletop.

Easy positioning

- Extra-wide arm board surface.
- Symmetrical design: may be positioned to support arm on either side of table.
- Includes pad and strap.