

100US

Gyn/Fertility Mobile Imaging Table
Aspect™ Series table line



Ideal C-Arm imaging solution for Hysterosalpingograms

- Carbon fiber tabletop with low x-ray attenuation
- Detachable leg support tabletop extension
- Unobstructed fluoroscopy imaging area: 22"W x 32"L (56 x 81 cm)
- Motorized tabletop elevation: 27" – 38" (66 – 96.5cm)
- 24"W x 46"L tabletop (61 x 117 cm)
- Patient weight capacity: 450 lbs. (204 kg)
- Standard US table side-rail: compatible with many lithotomy legholders and clamps
- 120 VAC operation with battery back-up operation (230VAC optional)
- Easy table transport with 5" locking casters
- Table weight = 375 lbs. (170 kg)
- 3-year parts warranty

Accessories included:

- Basic knee crutches with two side rail clamps
- Tabletop pad
- Foot control
- Patient restraint straps

Optional accessories available:

- Armboard
- Table hand control



Image Diagnostics, Inc.

310 Authority Drive Fitchburg, MA, 01420 USA Tel: 978-829-0009 Fax: 978-829-0027
www.imagediagnostics.com D000-517G

idi

Image Diagnostics, Inc.
Safety, Simplicity, Quality