

Image Diagnostics Launches Aspect ISR G3 Mobile Imaging Table

[SHARE](#) | [E-MAIL](#) | [PRINT](#) | [BOOKMARK](#) |

January 4, 2018—Image Diagnostics Inc. (IDI) announced the launch its Aspect ISR G3 imaging table. The float top carbon fiber table has 32 inches of longitudinal travel and up to 80 inches in artifact-free imaging area with a compact, efficient platform.

As a solution for C-arm imaging, the G3's four-way floating tabletop motion is operated via a table-mounted panning hand controller that operates all table functions including pan, up/down, Trendelenberg, lateral, and Isocentric roll. IDI's Isocentric lateral roll function allows the focused anatomy to be centered even when rolled laterally. The table also has no hydraulics or pneumatic drive systems, providing operational simplicity and reliability, noted the company.