

IDI Aspect 100 UCPlus Accessories & Consumables









idi

Image Diagnostics Inc.



Veteran
Owned





Aspect 100UCPlus Optional and Replacement Accessories:

Photo	Description Red text = Metric accessory rails	Catalog Number
	<p>IDI ComfortFin 350 Lithotomy Leg holders (Stirrups pair only) ComfortFin Lithotomy Leg holders are designed to improve infection control with durable, easy to clean silicone straps. Includes integrated lateral fin to help protect the head of the fibula and the peroneal nerve. 4" more working space with an increased range of motion. 350 lbs (159 kg).</p>	C100-3042
	<p>Kit: IDI ComfortFin350 Lithotomy Leg holders C100-3042 with C100-3046 economy side-rail clamps</p>	<p>A100-2246 A100-3074</p>
	<p>IDI ComfortFin 500 Lithotomy Leg holders (Stirrups pair only) ComfortFin Lithotomy Leg holders are designed to improve infection control with durable, easy to clean silicone straps. Includes integrated lateral fin to help protect the head of the fibula and the peroneal nerve. Boots are sized for today's heavier patients. 4" more working space with an increased range of motion. Patient weight capacity: 500 lbs. (227kg).</p>	C100-3045
	<p>Kit: IDI ComfortFin 500 Lithotomy Leg holders C100-3045 with C100-3046 economy side-rail clamps.</p>	<p>A100-2247 A100-3075</p>
	<p>Side-rail clamp, removable. IDI 4100. Attaches and locks anywhere on the side rail for fast and easy mounting . A push button safely releases rail clamp—helps to prevent false rail engagement. Autoclavable.</p>	<p>C100-3043 C100-3076</p>
	<p>Economy side-rail clamp. Attaches anywhere on the table rail for fast and easy mounting. Rounded ergonomic shape is easy on the hands. Autoclavable.</p>	<p>C100-3046 C100-3077</p>
	<p>Lithotomy Leg Holder Cart Easy to push/pull with ergonomic handles. Stirrups are isolated to prevent damage. Clamp storage compartment. Designed to safely store ComfortFin stirrups when not in use.</p>	C100-3044
	<p>Single Use Stirrup Boot Cover Waterproof single use stirrup boot covers help eliminate strike through and bacterial migration. Infectious patients (known or suspected). Patients with open sores, abrasion or lesions. Protect against cross contamination. <i>Available Q3 2019.</i></p>	C100-3070

Aspect 100UCPlus Optional and Replacement Accessories:

Photo	Description Red text = Metric accessory rails	Catalog Number
	Physician Elbow Supports (pair) Removable	C000-0594
	Armboard Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007 or A100-2379 above) unless used on a table with built-in accessory rails.	A310-056 A310-062
	Split Leg Holders For use with Aspect 100UC Plus Tables. Includes pads. Prone percutaneous, lateral or modified lateral nephrostomy. Includes two split leg holder clamps # C000-1115 or C000-1217 .	C000-1146 C000-1216
	Replacement Split Leg Holder Clamp For use with split leg holders. # C000-1146 and C000-1216	C000-1115 C000-1217
	IV Pole Stainless steel with four bottle hooks Range: 1070 -1660 mm Tube diameter: 25 / 18 mm IV Pole Clamp C000-1325 is required.	C000-1324
	IV Pole Clamp Aluminum. 25 mm pole opening	C000-1325
	Clamp-on Accessory Rail Rail is 8 ½" long x ¾" thick x 1-1/8" high. Clamp-on Accessory Rail 10mm x 25mm x 216mm clamp-on accessory rail.	A100-2218 A100-2838
	Stirrups Mounting Rail extension Kit Included with table. For use with leg-holders (Leg-holders are sold separately).	#A100-1876 A100-2676





Aspect 100UCPlus Optional and Replacement Accessories:

Photo	Description Red text = Metric accessory rails	Catalog Number
	Hand Control for 100UCPlus and 100UC. Controls height up/down, trendelenburg/revers, slide, cross table movements. (Included with table).	K000-0158
	Foot Switch Control for 100UCPlus. (Included with table).	A100-1540
	Rotating, with socket for attachment of accessories with 5/8 inch (15.88mm) circular post. Requires A100-1007 (or A100-2379 for metric rails) clamp-on accessory rail.	C000-0746 C000-1065
	Removable radiographic cassette tray. Loads from either side of table.	W100-1001

Aspect 100UCPlus Consumables

Photo	Description Red text = Metric accessory rails	Catalog Number
	Vinyl Urology Table Covers Box of 20	C000-0645
	Uro Drain Bags with Hose Disposable drain bag and hose assembly. Box of 20	Sterile: # C000-0593 Non-sterile: #C000-1111
	Disposable Covers for Aspect 100UCPlus Foot Controls 24" x 20" Box of 50	C000-0492
	Uro Drain Support Hoop Collapsible & removable drain bag support hoop.	C000-0595

Aspect 100UCPlus Consumables

Photo	Description Red text = Metric accessory rails	Catalog Number
	3"/ 76mm Restraint Straps (Pair) Fabric with velcro fasteners for Aspect 100UCPlus table.	C000-0571
	Restraint Straps One each 3" and 5" (76mm and 127mm) Fabric with velcro fasteners for Aspect 100UCPlus table.	C000-0328
	Uro Collection Unit 5 gallon/18.9 liter capacity capacity, Box of 10. Soft sided Vinyl container gently expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Built in zipper and handles allow for storage, closure of drain bag and hose assembly for enhanced infection control during transport and disposal.	C000-0612
	Uro Collection Unit 15 liter capacity, Box of 10. Flexible plastic sided container expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Integrated valve system allows user to shut off the flow of liquid from the drain bag in order to change collection units..	C000-1307

Aspect 100UCPlus Shipping Crate Specifications

- Weight (approx.): 815 lbs. (370kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

Image Diagnostics, Inc.

Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

Image Diagnostics, Inc.

310 Authority Drive
Fitchburg, Massachusetts
01420

For customer service:
sales@imagediagnostics.com

Toll free: 877.304.5434
Main: 978.829.0009
Fax: 978.829.0027

ImageDiagnostics.com

idi

Image Diagnostics Inc.



A Veteran-Owned Business