


ILEX C

Utility cart with accessory rail

The perfect compliment for motorized C-Arm control panels or rail-mounted accessory devices

Elegantly contoured, counter weighted, five caster (two lockable) dual-wheeled utility cart. Complete with adjustable U.S. or metric 12" (305MM) medical accessory rail and accessory basket.

- Two locking casters
- Utility basket with cable wraps
- 12" (305mm) adjustable medical accessory rail (U.S. or metric standard)
- Accessory rail is 42" (1067mm) from the floor
 - Control handle
 - Up to 30 lbs (13 kg) capacity
- Approximately 22" x 23" footprint



idi
Image Diagnostics, Inc.
Safety, Simplicity, Quality

Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging. Our product focus is in electro-mechanical solutions for fluoroscopic and digital and radiographic applications. IDI is committed to delivering superior product design, value and customer support in every product we build.

IDI IMAGING TABLES

IDI Cardiovascular Tables

IDI's Aspect line of interventional and endovascular tables deliver extended positioning and imaging capabilities.



IDI Urology Tables

For over 25 years IDI has manufactured the market leading mobile table for Urology.



IDI Pain Management Tables

C-arm tables for GI and pain management procedures. Tables with tilt, roll and longitudinal tabletop movements are available.



IDI VIDEO INTEGRATION SYSTEMS

ilex55 , ilex55T and ilex 58TC (Cardio) Video Integration Systems

The first mobile large field 4K UHD display solution with onboard video integration. The systems are equipped with 55" or 58" surgical or DICOM-ready high resolution UHD monitors. The T and TC Cardio systems feature a touch-screen video management system providing image layouts similar to those in hybrid operating rooms and state-of-the-art Cath labs.



MDS Video Integration System

Advanced Video Integration for multi-modality image visualization in Urology, GI & Vascular surgery.



ilex32 Monitor Stand

Portable image display when you just need a bigger monitor.



idi

Image Diagnostics, Inc.
Safety, Simplicity, Quality



Veteran-owned

*Please reference IDI Warranty Document D000-2069
Image Diagnostics, Inc. 310 Authority Drive Fitchburg, MA, 01420 USA
Tel: 978-829-0009 Fax: 978-829-0027 imagediagnostics.com L280-0940 Rev C