

# Connexus Direct Connect Q-CARD

**idi**  
Image Diagnostics, Inc.



## **WARNING!**

**Read operator manual prior to use**



### **Video Panel**

- 1) Each channel has two connections DVI or SOG (Sync on Green)
- 2) Only one video source can be transmitted at a time, per channel
- 3) Route cables between Connexus and cable relief bar

### **Cleaning and Storage:**

- 1) Unplug power cord and move to safe location
- 2) Wipe Down Surface
- 3) Clean and disinfect, see manual for approved cleaners

### **AC Power:**

- 1) Plug power cord into properly grounded outlet
- 2) Protect cord from damage by casters or other equipment



Scan to  
Access  
User  
Manuals  
Online