



ASPECT

# 100US

*Gyn/Fertility Mobile Imaging Table*



### **Ideal C-Arm imaging solution for Hysterosalpingograms**

- Unobstructed imaging area: 21"W x 32"L (53 x 81 cm)
  - Carbon fiber tabletop - 24"W x 46"L tabletop (61 x 117 cm)
  - Detachable boarding extension – extremity support only
  - Motorized 27" – 38" (66 – 96.5cm)
  - Patient weight capacity: 450 lbs. (204 kg)
  - Standard US side-rails:
    - Compatible with many lithotomy leg holders and clamps
  - 120 VAC operation with battery back-up operation
  - Large 5" locking casters
  - Table weight = 375 lbs. (170 kg)
- Warranty: One year labor/ three years parts\*

### **Accessories included with package A100-1479**

- Basic knee crutches with two side rail clamps
- Tabletop pad
- Foot control
- Patient restraint straps

### **Optional accessories available:**

- Armboard
- Table hand control



**idi**  
imagediagnostics.com

**INTELLIGENT DESIGN ► TANGIBLE OUTCOMES**

# IDI IMAGING TABLES

## IDI Interventional / Endovascular Tables

IDI's Aspect line of interventional and endovascular tables deliver extended positioning and imaging capabilities.



## IDI Urology / GYN Tables

For over 25 years IDI has manufactured the market leading mobile table for Urology.



## IDI Spinal Imaging Pain Management Tables

C-arm tables for GI and pain management procedures. Tables with tilt, roll and longitudinal tabletop movements are available.



# IDI VIDEO INTEGRATION SYSTEMS

## Ilex Video Integration Systems

The first mobile large field 4K UHD display solution with built in integration. A State-of-the-Art Display Solution for Vascular, GI and Neurosurgery.



## MDS Video Integration System

Advanced Video Integration for multi-modality image visualization in Urology, GI & Vascular surgery.



## Monitor Carts

Portable image displays when you need a bigger monitor or a bigger monitor with recording capabilities. Available in both HD and UHD 4K versions.



Image Diagnostics, Inc. 310 Authority Drive Fitchburg, MA 01420 Tel: 978-829-0009  
A veteran-owned company since 1986. Designed and manufactured in the U.S.  
Image Diagnostics, Inc. reserves the right to update product specifications without notification.  
\*Warranty information for IDI products is available upon request from Image Diagnostics, Inc.  
L100-3115 Rev B

★★★★★  
Veteran-Owned  
Assembled in  
the USA with  
global materials  
★★★★★

