

Image Diagnostics, Inc.

BOM NUMBER :A240-1140 REV: C
BOM DESCRIPTION:KIT, HD MONITOR MOUNTING HARDWARE
PRINT DATE :05/11/2022 09:07:58 AM
EFFECTIVE DATE :05/11/2022
QTY :1.0000

BOM STATUS: Released
MODEL NO:

PAGE NO:1
UNCSIN.FRX

UNCOSTED SINGLE LEVEL BOM REPORT

ITEM	BOM/PART NUMBER	REV	DESCRIPTION	MANUFACTURER	MODEL NUMBER	QTY PER ASSY	QTY REQUIRED	COST	COMMENTS
0001	C000-0931	A	SPACER,SHDR,.125 LG,NYL,5/16 SC,3/8	MCMaster	91145A178	4.0000	4.0000		
0002	C000-1484	A	SPACER,SHDR,NYL,.26IDx.38-.56ODx.12	McMaster CARR	90062A015	6.0000	6.0000		
0003	C000-1485	A	SCR,SHCS,ALY,6MM X 18MM			6.0000	6.0000		
0004	C000-1486	A	WSHR,FLT,ALY,6MM			6.0000	6.0000		
0005	C000-1487	A	WSHR,LCK,ALY,1/4IN			6.0000	6.0000		
0006	C000-1499	A	SCR,SHCS,ALY,8MM X 20MM			4.0000	4.0000		
0007	C000-1500	A	WSHR,FLT,ZP,8MM			4.0000	4.0000		
0008	C000-1501	A	WSHR,LCK,ZP,8MM			4.0000	4.0000		
0009	Z240-1138	A	TAPE,POLYETH.040X3/8X5 CUT TO SIZE			8.0000	8.0000		