# IDI 32" 4K MONITOR CART

# Designed for use as a large table side monitor

- Direct connect to C-arm workstations and cameras
- Digital connectivity with no scaling or conversion required\*
- · Optional HD or UHD recording

The IDI Monitor Cart provides a cost effective solution for those customers requiring an additional monitor with mirrored images in multiple locations.

### Cart includes:

- Single or split-screen image display
- Heavy duty cart stand
- Maximum height: 72.5"/184cm
- ±10° monitor tilt
- · Cable channel
- Hygienic non-porous work surface with nook
- 15 foot power cord (NEMA 5-15PH)
- · Cable hook
- Cart includes video input cables
- One year parts and labor warranty
- Lower shelf designed to enclose and support a Sony Printer (Optional)

#### **Monitor**

- 32"/800mm 4K medical grade monitor with medical bios
- 3840x2160 UHD 4K
- 1350:1 Contrast ratio
- \*3G/12G SDI upgrade option available

## Input options

- HDMI, 2 x DisplayPort, DVI
- Optional SOG is required to connect to legacy C-Arm systems



Monitor Cart with 32" 4K display

\*Standard inputs and frame rates Monitor specifications are subject to change, please contact Image Diagnostics for latest monitor specifications



ORDER: Sales@ImageDiagnostics.com

SUPPORT: TechSupport@ImageDiagnostics.com

WEBSITE: ImageDiagnostics.com

+1-978-829-0009 Monday through Friday 8am to 5pm, Eastern



