

Interventional / Endovascular Imaging

ISO-DRIVE

Table and Accessories Catalog



idi
imagediagnostics.com

INTELLIGENT DESIGN ► TANGIBLE OUTCOMES

Iso-Drive Advanced Spinal Image Table Options and Accessories Catalog



Iso-Drive Description:

C-Arm table for pain management and other imaging procedures where motorized elevation, tilt, lateral roll and longitudinal motion are required.

- Isocentric Roll
- Elevation
- Trendelenburg tilt

Iso-Drive Spinal Imaging Table #A100-3121 with .375" x 1.125" accessory rails (120VAC - 240VAC)

Includes the following:

- X100-3849 Tabletop Pad, 2" (51mm) thick
- X100-3212: Facial Cushion, 1" Thick
- X100-3234: Head Cushion, 3.375" Thick
- X100-3235: Facial Cushion, 3.375" Thick
- X100-3236: Tapered Facial Cushion, 3.375" Thick
- K000-0299: Hand control
- A100-3731: Patient restraint straps (2 straps, 1 pad)
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 24" (610mm) x 84" (2134mm)
with facial cutout







Iso-Drive General Purpose Table #A100-3120 (120 - 240 VAC)

Includes the following:







- X100-3310: Tabletop pad, 2" (51mm) thick
- K000-0299: Hand control
- A100-3731: Patient restraint straps (2 straps, 1 pad)
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 24" (610mm) x 84" (2134mm)

Iso-Drive Table Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Replacement Tabletop Pad For A100-3120 table. 2" thick foam pad 24" wide x 84" long.	X100-3310
	Replacement Tabletop Pad For A100-3121 table. 2" thick for 84" x 24" x 19" top with 12" facial relief and cutout.	X100-3849
	Snap-on Socket Accommodates a round accessory post. Quickly snaps onto rail and pull bottom lever to release. Mounts and locks anywhere on standard OR table rail. Lithotomy wheel for accessory angle adjustments.	C000-0746 C000-1065
	Simple clamp Accommodates a round accessory post. All metal design clamp mounts and locks anywhere on standard OR table rail. Autoclavable.	C100-3046
	Snap-on Blade Clamp Accommodates a blade-style accessory post. Quickly snaps onto rail, push button to release clamp. Mounts and locks anywhere on standard OR table rail. Prevents false rail engagement. Autoclavable.	C100-3043
	Simple Blade Clamp Accommodates blade-style accessory post. Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost-effective solution.	C000-1115

Iso-Drive Table Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Armboard Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007 or A100-2379 .	A310-056 A310-062
	Replacement pad for A310-056 or A310-062 Armboard	C100-3288
	Single use soft armband straps, package of 12	C000-1555
	Clamp-on Accessory Rail Rail is 8.5" long x .375" thick x 1.25" high. Clamp-on Accessory Rail 10mm x 25mm x 216mm clamp-on accessory rail.	A100-1007 A100-2379
	Arm Board Carbon fiber, shoulder mount (attaches by sliding thin section under patient / table pad), flat. Includes pad and strap.	A310-059
	Replacement Pad for A310-059 Armboard	X100-3244
	Hand Control, Iso-Drive Table Must verify table serial # with IDI for compatibility. Motions controlled: Elevation, trendelenburg, Iso-Roll. (Included with table).	K000-0299
	Scatter-Guard radiation shield draping - .5mm lead equivalence For 24" (610mm) rectangular tabletop X100-995. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop.	A610-096
	Radiation Shield - .5mm lead equivalence Fixed shield with integrated clamp system. Clamps to tables that do not have side rails. Snap-on upper section above level of tabletop is removable. 24" wide x 32.125" high x 4" deep.	A610-080
	<ul style="list-style-type: none"> • Facial Cushion, crescent shape. • Disposable covers for facial cushion, box of 50. 	Cushion: #C000-0597 Covers: #C000-0598

Iso-Drive Table Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Heavy Duty Medical Grade Anti-Fatigue Mat (24" (610mm) x 66" (1676mm) x .75" (19mm) IDI heavy-duty anti-fatigue mat. One-piece construction for high-traffic areas. Safe & non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.	C000-1181
	IV Pole (C000-1324) Stainless steel with four bottle hooks and IV Pole Clamp (C000-1325) Range: 1070 -1660 mm Tube diameter: 25 / 18 mm	A100-3285
	Compass Contour Drape Custom shaped drape for the Compass Contour Radiation Shield. Slips over the lower portion of the shield and is held in place with an elastic band. Sterile. Case of 20.	X610-0208
	Curve Spine Frame Provides maximum lordosis for Laminectomy, Decompression, Disc Surgery and Microdiscectomy procedures. Includes hand crank. Patient supports are adjustable to fit different patient sizes.	C000-1750
	Compass Contour Radiation Shield Contoured ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~30" (76cm) x 23.5" (60cm) with a 10" radius along the bottom edge of the leaded acrylic shield. Does not include ceiling suspension arm.	A610-0200
	Compass Frame Radiation Shield Rectangular ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~ 30" (76cm) x 23.5" (60cm). Does not include ceiling suspension arm.	A610-0215

Iso-Drive Table Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Restraint Straps (Pair) Set of two patient restraint straps with buckles and one pad.	A100-3731
ISO-DRIVE Shipping Crate Specifications		
Weight (approx.): 620 Lb (281 kg.) Length: 90" (2,286mm) Width: 32" (813mm) Height: 32" (813mm)		



Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

Image Diagnostics, Inc.
310 Authority Drive
Fitchburg, Massachusetts
01420

For customer service:
sales@imagediagnostics.com
Telephone: 978.829.0009

INTELLIGENT DESIGN ► TANGIBLE OUTCOMES



ORDER: sales@imagediagnostics.com
SUPPORT: techsupport@imagediagnostics.com
WEBSITE: <https://imagediagnostics.com>

+1-978-829-0009
Monday through Friday

