



100US

Gyn/Fertility Mobile Imaging Table



Ideal C-Arm imaging solution for Hysterosalpingograms

- Unobstructed imaging area: 21"W x 32"L (53 x 81 cm)
 - Carbon fiber tabletop - 24"W x 46"L tabletop (61 x 117 cm)
 - Detachable boarding extension – extremity support only
 - Motorized 27" – 38" (66 – 96.5cm)
 - Patient weight capacity: 450 lbs. (204 kg)
 - Standard US side-rails:
 - Compatible with many lithotomy leg holders and clamps
 - 120 VAC operation with battery back-up operation
 - Large 5" locking casters
 - Table weight = 375 lbs. (170 kg)
- Warranty: One year labor/ three years parts*

Accessories included with package A100-1479

- Basic knee crutches with two side rail clamps
- Tabletop pad
- Foot control
- Patient restraint straps

Optional accessories available:

- Armboard
- Table hand control



idi
imagediagnostics.com

INTELLIGENT DESIGN ► TANGIBLE OUTCOMES

IDI IMAGING TABLES

IDI Interventional / Endovascular Tables

IDI's Aspect line of interventional and endovascular tables deliver extended positioning and imaging capabilities.



IDI Urology / GYN Tables

For over 25 years IDI has manufactured the market leading mobile table for Urology.



IDI Spinal Imaging Pain Management Tables

C-arm tables for GI and pain management procedures. Tables with tilt, roll and longitudinal tabletop movements are available.



IDI VIDEO INTEGRATION SYSTEMS

Ilex Video Integration Systems

The first mobile large field 4K UHD display solution with built in integration. A State-of-the-Art Display Solution for Vascular, GI and Neurosurgery.



MDS Video Integration System

Advanced Video Integration for multi-modality image visualization in Urology, GI & Vascular surgery.



Monitor Carts

Portable image displays when you need a bigger monitor or a bigger monitor with recording capabilities. Available in both HD and UHD 4K versions.



Image Diagnostics, Inc. 310 Authority Drive Fitchburg, MA 01420 Tel: 978-829-0009
A veteran-owned company since 1986. Designed and manufactured in the U.S.
Image Diagnostics, Inc. reserves the right to update product specifications without notification.
*Warranty information for IDI products is available upon request from Image Diagnostics, Inc.
L100-3115 Rev B

★★★★★
Veteran-Owned
Assembled in
the USA with
global materials
★★★★★

