

## General Imaging / Pain Management

# 100RTL 100RT 100T 100E Tables and Accessories



**100RTL** Motorized elevation, longitudinal travel, lateral and Trendelenburg tilt



**100RT** Motorized elevation, lateral and Trendelenburg tilt



**100T** Motorized elevation and Trendelenburg tilt



**100E** Motorized elevation



## 100RTL Table



### 100RTL Description:


C-Arm table for pain management and other imaging procedures where motorized elevation, tilt, lateral roll and longitudinal motion are required.

- Motorized elevation
- Trendelenburg tilt
- Lateral roll
- Longitudinal motion
- 120VAC is standard
- 230VAC is optional and must be specified at time of order

### 100RTL #A100-1254 (120VAC)    A100-2637 (230VAC)

Includes the following:

- X100-1742: Tabletop pad, 2" (51mm) thick
- K000-0197: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial




Tabletop 24" (610mm) x 80" (2032mm)

### 100RTL #A100-2646 (120VAC)    A100-2647 (230VAC)

Includes the following:

- X100-1757: Tabletop pad, 2" (51mm) thick
- K000-0197: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 22" (559mm) x 80" (2032mm)

## 100RTL #A100-2672 (120VAC) # A100-2673 (230VAC)








Includes the following:

- X100-2083: Tabletop pad, 2" (51mm) thick
- K000-0197: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion




Tabletop 24" (610mm) x 80" (2032mm)









### 100RTL Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<b>Replacement Tabletop Pad</b> For A100-1254 or A100-2637 tables. 2" thick foam pad 24" wide x 80" long.	X100-1742
	<b>Replacement Tabletop Pad</b> For A100-2646 or A100-2647 tables. 2" thick foam pad. 22" wide with a facial cutout x 80" long.	X100-1757
	<b>Replacement Tabletop Pad</b> For A100-2672 or A100-2673 tables. 2" thick foam pad 28" wide x 80" long.	X100-2083
	<b>Armboard</b> Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007 or A100-2379).	A310-056 <span style="color: red;">A310-062</span>
	<b>Replacement pad for A310-056 or A310-062 Armboard</b>	C100-3288
	<b>Single use soft armband straps, package of 12</b>	C000-1555
	<b>Clamp-on Accessory Rail</b> Rail is 8.5" long x .375" thick x 1.25" high. <span style="color: red;">Clamp-on Accessory Rail</span> <span style="color: red;">10mm x 25mm x 216mm clamp-on accessory rail.</span>	A100-1007 <span style="color: red;">A100-2379</span>
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient / table pad), flat. Includes pad and strap.	A310-059
	<b>Replacement Pad for A310-059 Armboard</b>	X100-3244
	<b>Hand Control, 100RTL</b> (must verify table serial # with IDI for compatibility). Motions controlled: Elevation, trendelenburg, roll, slide. (Included with table).	K000-0197




## 100RTL Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Foot Control for 100RTL</b> Foot Control unit controls all table motions.</p>	A100-1961
	<p><b>Scatter Guard radiation shield draping - .5mm lead equivalence</b> For 28" (711mm) wide optional tabletop X100-1325 Dimensions: 59" (1500mm) long, 24" (610mm) below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop.</p>	A610-108
	<p><b>Radiation Shield - .5mm lead equivalence</b> Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables. 24" wide x 32.125" high x 4" deep.</p>	A610-080
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b> For for 24" (610mm) rectangular tabletop X100-995. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop.</p>	A610-096
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b> for 22" (559mm) tabletop X100-733 models with 22" (559mm) tapered-end tabletop with facial cutout. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop.</p>	A610-092
	<ul style="list-style-type: none"> <li>• <b>Facial Cushion, crescent shape.</b></li> <li>• <b>Disposable covers for facial cushion, box of 50.</b></li> </ul>	<p>Cushion: #C000-0597 Covers: #C000-0598</p>

## 100RTL Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Curve Spine Frame</b> Provides maximum lordosis for Laminectomy, Decompression, Disc Surgery and Microdiscectomy procedures. Includes hand crank. Patient supports are adjustable to fit different patient sizes.</p>	C000-1750
	<p><b>Heavy Duty Medical Grade Anti-Fatigue Mat</b> (24" (610mm) x 66" (1676mm) x .75" (19mm) IDI heavy-duty anti-fatigue mat. One-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181
	<p><b>Compass Contour Radiation Shield</b> Contoured ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~30" (76cm) x 23.5" (60cm) with a 10" radius along the bottom edge of the leaded acrylic shield. Does not include ceiling suspension arm.</p>	A610-0200
	<p><b>Compass Frame Radiation Shield</b> Rectangular ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~ 30" (76cm) x 23.5" (60cm). Does not include ceiling suspension arm.</p>	A610-0215
	<p><b>Compass Contour and Frame Drape</b> Custom shaped drape for the Compass Contour and Frame Radiation Shield. Slips over the entire shield and is held in place with an elastic band. Sterile. Case of 20.</p>	C610-0256
	<p><b>IV Pole C000-1324</b> Stainless steel with four bottle hooks Range: 1070 -1660 mm Tube diameter: 25 / 18 mm <b>IV Pole Accessory Rail Clamp C000-1325 included</b></p>	A100-3285
	<p><b>Simple clamp</b> Accommodates a round accessory post. All metal design clamp mounts and locks anywhere on standard OR table rail. Autoclavable.</p>	C100-3046
	<p><b>Snap-on Blade Clamp</b> Accommodates a blade-style accessory post. Quickly snaps onto rail, push button to release clamp. Mounts and locks anywhere on standard OR table rail. Prevents false rail engagement. Autoclavable.</p>	C100-3043

## 100RTL Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Simple Blade Clamp</b> Accommodates blade-style accessory post. Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost-effective solution.</p>	C000-1115
	<p><b>Snap-on Socket</b> Accommodates a round accessory post. Quickly snaps onto rail and pull bottom lever to release. Mounts and locks anywhere on standard OR table rail. Lithotomy wheel for accessory angle adjustments.</p>	C000-0746
	<p><b>Restraint Straps</b> Set of two patient restraint straps with buckle locking mechanism and one pad for patient safety and comfort</p>	A100-3731

## 100RTL Shipping Crate Specifications

- Weight (approx.): 605 Lb. (274kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

# 100RT Mobile Imaging Table



## 100RT Description:

C-arm table for pain management and other imaging procedures where motorized elevation, tilt and lateral roll are required.

- Motorized Elevation
- Trendelenburg tilt
- Lateral Roll
- 120VAC is standard
- 230VAC is optional and must be specified at time of order

## 100RT #A100-870 (120VAC)    A100-2636 (230VAC)

Includes the following:

- X100-1742: Tabletop pad, 2" (51mm) thick
- K000-0195: Hand control.
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 24" (610mm) x 80" (2032mm)

## 100RT #A100-2644 (120VAC)    A100-2645 (230VAC)

Includes the following:

- X100-1757: Tabletop pad, 2" (51mm) thick
- K000-0195: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 22" (610mm) x 80" (2032mm)








**100RT #A100-2660 (120VAC) #A100-2661(230VAC)**

Includes the following:

- X100-2083: Tabletop pad, 2" (51mm) thick
- K000-0195: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50

Tabletop 28" (711mm) x 80" (2032mm)

**100RT Table Optional and Replacement Accessories:**







Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<b>Replacement Tabletop Pad</b> For A100-870 or A100-2636 tables. 2" thick foam pad 24" wide x 80" long.	X100-1742
	<b>Replacement Tabletop Pad</b> For A100-2644 or A100-2645 tables. 2" thick foam pad. 22" wide with a facial cutout x 80" long.	X100-1757
	<b>Replacement Tabletop Pad</b> For A100-2660 or A100-2661 tables. 2" thick foam pad 28" wide x 80" long.	X100-2083
	<b>Armboard</b> Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007 or A100-2379).	A310-056 <span style="color: red;">A310-062</span>
	<b>Replacement pad for A310-056 or A310-062 Armboard</b>	X100-1755
	<b>Clamp-on Accessory Rail</b> Rail is 8.5" long x .375" thick x 1.25" high. <span style="color: red;">Clamp-on Accessory Rail</span> <span style="color: red;">10mm x 25mm x 216mm clamp-on accessory rail.</span>	A100-1007 <span style="color: red;">A100-2379</span>
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient / table pad), flat. Includes pad and strap.	A310-059
	<b>Replacement Pad for A310-059 Armboard</b>	X100-3244
	<b>Hand Control, 100RT</b> Must verify table serial # with IDI for compatibility. (Included with table).	K000-0195



## 100RT Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Foot Control</b> For 100RT table. Must verify table serial # with IDI for compatibility.</p>	A100-1960
	<p><b>Scatter Guard radiation shield draping - .5mm lead equivalence</b> For 28" (711mm) wide optional tabletop X100-1325 Dimensions: 59" (1500mm) long, 24" (610mm) below tabletop. Draping extends along both sides and head-end of table. Supported by fabric section over tabletop.</p>	A610-108
	<p><b>Radiation Shield - .5mm lead equivalence</b> Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables. 24" wide x 32.125" high x 4" deep.</p>	A610-080
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b> For 24" (610mm) rectangular tabletop X100-995. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop</p>	A610-096
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b> for 22" (559mm) tabletop X100-733 models with 22" (559mm) tapered-end tabletop with facial cutout. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop.</p>	A610-092
	<ul style="list-style-type: none"> <li>• <b>Facial Cushion, crescent shape.</b></li> <li>• <b>Disposable covers for facial cushion, box of 50.</b></li> </ul>	<p>Cushion: #C000-0597 Covers: #C000-0598</p>

## 100RT Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Curve Spine Frame</b> Provides maximum lordosis for Laminectomy, Decompression, Disc Surgery and Microdiscectomy procedures. Includes hand crank. Patient supports are adjustable to fit different patient sizes.</p>	C000-1750
	<p><b>Heavy Duty Medical Grade Anti-Fatigue Mat</b> (24" (610mm) x 66" (1676mm) x .75" (19mm) IDI heavy-duty anti-fatigue mat. One-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181
	<p><b>Compass Contour Radiation Shield</b> Contoured ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~30" (76cm) x 23.5" (60cm) with a 10" radius along the bottom edge of the leaded acrylic shield. Does not include ceiling suspension arm.</p>	A610-0200
	<p><b>Compass Frame Radiation Shield</b> Rectangular ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~ 30" (76cm) x 23.5" (60cm). Does not include ceiling suspension arm.</p>	A610-0215
	<p><b>Compass Contour and Frame Drape</b> Custom shaped drape for the Compass Contour and Frame Radiation Shield. Slips over the entire shield and is held in place with an elastic band. Sterile. Case of 20.</p>	C610-0256
	<p><b>IV Pole C000-1324</b> Stainless steel with four bottle hooks Range: 1070 -1660 mm Tube diameter: 25 / 18 mm <b>IV Pole Accessory Rail Clamp C000-1325 included</b></p>	A100-3285

## 100RT Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Simple clamp</b> Accommodates a round accessory post. All metal design clamp mounts and locks anywhere on standard OR table rail. Autoclavable.</p>	C100-3046
	<p><b>Snap-on Blade Clamp</b> Accommodates a blade-style accessory post. Quickly snaps onto rail, push button to release clamp. Mounts and locks anywhere on standard OR table rail. Prevents false rail engagement. Autoclavable.</p>	C100-3043
	<p><b>Simple Blade Clamp</b> Accommodates blade-style accessory post. Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost-effective solution.</p>	C000-1115
	<p><b>Snap-on Socket</b> Accommodates a round accessory post. Quickly snaps onto rail and pull bottom lever to release. Mounts and locks anywhere on standard OR table rail. Lithotomy wheel for accessory angle adjustments.</p>	C000-0746
	<p><b>Restraint Straps</b> Set of two patient restraint straps with buckle locking mechanism and one pad for patient safety and comfort</p>	A100-3731

## 100RT Shipping Crate Specifications

- Weight (approx.): 571 Lb(259 kg.)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

# 100T Mobile Imaging Table



## 100T Description:


C-arm table for pain management and other procedures where motorized elevation and Trendelenburg tilt are required.

- Motorized Elevation
- Trendelenburg tilt
- 120VAC is standard
- 230VAC is optional and must be specified at time of order

## 100T #A100-658 (120VAC)    A100-2638 (230VAC)

Includes the following:

- X100-1742: Tabletop pad, 2" (51mm) thick
- K000-0054: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50




Tabletop 24" (610mm) x 80" (2032mm)

## 100T #A100-2648 (120VAC)    A100-2649 (230VAC)

Includes the following:

- X100-1757: Tabletop pad, 2" (51mm) thick
- K000-0054: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50




Tabletop 22" (559mm) x 80" (2032mm)








## 100T #A100-2658 (120VAC)    A100-2659 (230VAC)

Includes the following:

- X100-2083: Tabletop pad, 2" (51mm) thick
- K000-0054: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50









### 100T Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<b>Replacement Tabletop Pad</b> For A100-658 or A100-2638 tables. 2" thick foam pad 24" wide x 80" long.	X100-1742
	<b>Replacement Tabletop Pad</b> For A100-2648 or A100-2649 tables. 2" thick foam pad. 22" wide with a facial cutout x 80" long.	X100-1757
	<b>Replacement Tabletop Pad</b> For A100-2658 or A100-2659 tables. 2" thick foam pad 28" wide x 80" long.	X100-2083
	<b>Armboard</b> Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007 or A100-2379.	A310-056 <span style="color: red;">A310-062</span>
	<b>Replacement pad for A310-056 or A310-062 Armboard</b>	C100-3288
	<b>Single use soft armband straps, package of 12</b>	C000-1555
	<b>Clamp-on Accessory Rail</b> Rail is 8.5" long x .375" thick x 1.25" high. <b>Clamp-on Accessory Rail</b> <span style="color: red;">10mm x 25mm x 216mm clamp-on accessory rail.</span>	A100-1007 <span style="color: red;">A100-2379</span>
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient / table pad), flat. Includes pad and strap.	A310-059
	<b>Replacement Pad for A310-059 Armboard</b>	X100-3244
	<b>Hand Control</b> for 100T table. Elevate and tilt functions. (Included with table).	K000-0054

## 100T Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Foot Control for 100T</b> Substitute for hand control. Cannot use hand and foot control at the same time.</p>	A100-500
	<p><b>Scatter Guard radiation shield draping - .5mm lead equivalence</b> For 28" (711mm) wide optional tabletop X100-1325 Dimensions: 59" (1500mm) long, 24" (610mm) below tabletop. Draping extends along both sides and head--end of table. Supported by fabric section over tabletop.</p>	A610-108
	<p><b>Radiation Shield - .5mm lead equivalence</b> Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables. 24" wide x 32.125" high x 4" deep.</p>	A610-080
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b> For 24" (610mm) rectangular tabletop X100-995. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop</p>	A610-096
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b> for 22" (559mm) tabletop X100-733 models with 22" (559mm) tapered-end tabletop with facial cutout. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop.</p>	A610-092
	<ul style="list-style-type: none"> <li>• <b>Facial Cushion, crescent shape.</b></li> <li>• <b>Disposable covers for facial cushion, box of 50.</b></li> </ul>	<p>Cushion: #C000-0597 Covers: #C000-0598</p>

## 100T Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Curve Spine Frame</b> Provides maximum lordosis for Laminectomy, Decompression, Disc Surgery and Microdiscectomy procedures. Includes hand crank. Patient supports are adjustable to fit different patient sizes.</p>	C000-1750
	<p><b>Heavy Duty Medical Grade Anti-Fatigue Mat</b> (24" (610mm) x 66" (1676mm) x .75" (19mm) IDI heavy-duty anti-fatigue mat. One-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181
	<p><b>Compass Contour Radiation Shield</b> Contoured ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~30" (76cm) x 23.5" (60cm) with a 10" radius along the bottom edge of the leaded acrylic shield. Does not include ceiling suspension arm.</p>	A610-0200
	<p><b>Compass Frame Radiation Shield</b> Rectangular ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~ 30" (76cm) x 23.5" (60cm). Does not include ceiling suspension arm.</p>	A610-0215
	<p><b>Compass Contour and Frame Drape</b> Custom shaped drape for the Compass Contour and Frame Radiation Shield. Slips over the entire shield and is held in place with an elastic band. Sterile. Case of 20.</p>	C610-0256
	<p><b>IV Pole C000-1324</b> Stainless steel with four bottle hooks Range: 1070 -1660 mm Tube diameter: 25 / 18 mm <b>IV Pole Accessory Rail Clamp C000-1325 included</b></p>	A100-3285

## 100T Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Simple clamp</b> Accommodates a round accessory post. All metal design clamp mounts and locks anywhere on standard OR table rail. Autoclavable.</p>	C100-3046
	<p><b>Snap-on Blade Clamp</b> Accommodates a blade-style accessory post. Quickly snaps onto rail, push button to release clamp. Mounts and locks anywhere on standard OR table rail. Prevents false rail engagement. Autoclavable.</p>	C100-3043
	<p><b>Simple Blade Clamp</b> Accommodates blade-style accessory post. Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost-effective solution.</p>	C000-1115
	<p><b>Snap-on Socket</b> Accommodates a round accessory post. Quickly snaps onto rail and pull bottom lever to release. Mounts and locks anywhere on standard OR table rail. Lithotomy wheel for accessory angle adjustments.</p>	C000-0746
	<p><b>Restraint Straps</b> Set of two patient restraint straps with buckle locking mechanism and one pad for patient safety and comfort</p>	A100-3731

## 100T Shipping Crate Specifications

- Weight (approx): 540 Lb /254kg
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)



# 100E Mobile Imaging Table



## 100E Description:

C-arm table for pain management and other procedures where motorized elevation is the only movement required.

- 120VAC is standard
- 230VAC is optional and must be specified at time of order

## 100E #A100-402 (120VAC)    A100-2635 (230VAC)

Includes the following:

- X100-1742: Tabletop pad, 2" (51mm) thick
- K000-0041: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion.
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 24" (610mm) x 80" (2032mm)

## 100E #A100-2642 (120VAC)    A100-2643 (230VAC)

Includes the following:

- X100-1757: Tabletop pad, 2" (51mm) thick
- K000-0041: hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 22" (559mm) x 80" (2032mm)

**100E #A100-2656 (120VAC)    A100-2657 (230VAC)**

Includes the following:

- X100-2083: Tabletop pad 2" (51mm) thick
- K000-0041: Hand control
- A100-3731: Patient restraint straps
- C000-0597: Crescent shaped facial cushion
- C000-0598: Disposable covers for crescent-shaped facial cushion, box of 50



Tabletop 28" (711mm) x 80" (2032mm)







**100E Table Optional and Replacement Accessories:**

Photo	Description <b>Red text = Metric accessory rails</b>	Catalog Number
	<b>Replacement Tabletop Pad</b> For A100-402 or A100-2635 tables. 2" thick foam pad 24" wide x 80" long.	X100-1742
	<b>Replacement Tabletop Pad</b> For A100-2642 or A100-2643 tables. 2" thick foam pad. 22" wide with a facial cutout x 80" long.	X100-1757
	<b>Replacement Tabletop Pad</b> For A100-2656 or A100-2657 tables. 2" thick foam pad 28" wide x 80" long.	X100-2083
	<b>Armboard</b> Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007 or A100-2379).	A310-056 <b>A310-062</b>
	<b>Replacement pad for A310-056 or A310-062 Armboard</b>	C100-3288
	<b>Single use soft armband straps, package of 12</b>	C000-1555
	<b>Clamp-on Accessory Rail</b> Rail is 8.5" long x .375" thick x 1.25" high. <b>Clamp-on Accessory Rail</b> 10mm x 25mm x 216mm clamp-on accessory rail.	A100-1007 <b>A100-2379</b>
	<b>Arm Board</b> Carbon fiber, shoulder mount (attaches by sliding thin section under patient / table pad), flat. Includes pad and strap.	A310-059
	<b>Replacement Pad for A310-059 Armboard</b>	X100-3244
	<b>Hand control for 100E.</b> (Included with table).	K000-0041

## 100E Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Foot Control 100E</b></p>	<p>A100-758</p>
	<p><b>Scatter Guard radiation shield draping - .5mm lead equivalence</b>            For 28" (711mm) wide optional tabletop X100-1325            Dimensions: 59" (1500mm) long, 24" (610mm) below tabletop. Draping extends along both sides and head-end of table. Supported by fabric section over tabletop.</p>	<p>A610-108</p>
	<p><b>Radiation Shield - .5mm lead equivalence</b>            Fixed shield with integrated clamp system. For use with tables that do not have side rails. Clamps to table sides or at head-end of table. Snap-on upper section above level of tabletop is removable. Use with 100RTL, 100RT, 100T, 100E tables. 24" wide x 32.125" high x 4" deep.</p>	<p>A610-080</p>
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b>            For 24" (610mm) rectangular tabletop X100-995.            Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop.</p>	<p>A610-096</p>
	<p><b>Scatter-Guard radiation shield draping - .5mm lead equivalence</b>            for 22" (559mm) tabletop X100-733 models with 22" (559mm) tapered-end tabletop with facial cutout. Dimensions: 59" (1500mm) long. 24" (610mm) draping below tabletop. Draping extends along both side and head-end of table. Supported by fabric section over tabletop.</p>	<p>A610-092</p>
	<ul style="list-style-type: none"> <li>• <b>Facial Cushion, crescent shape.</b></li> <li>• <b>Disposable covers for facial cushion, box of 50.</b></li> </ul>	<p>Cushion: #C000-0597 Covers: #C000-0598</p>

## 100E Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Curve Spine Frame</b> Provides maximum lordosis for Laminectomy, Decompression, Disc Surgery and Microdiscectomy procedures. Includes hand crank. Patient supports are adjustable to fit different patient sizes.</p>	C000-1750
	<p><b>Heavy Duty Medical Grade Anti-Fatigue Mat</b> (24" (610mm) x 66" (1676mm) x .75" (19mm) IDI heavy-duty anti-fatigue mat. One-piece construction for high-traffic areas. Safe &amp; non-toxic, contains no plasticizers and will not leach over time. Durable brushed texture.</p>	C000-1181
	<p><b>Compass Contour Radiation Shield</b> Contoured ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~30" (76cm) x 23.5" (60cm) with a 10" radius along the bottom edge of the leaded acrylic shield. Does not include ceiling suspension arm.</p>	A610-0200
	<p><b>Compass Frame Radiation Shield</b> Rectangular ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhanced eye, head and neck protection from x-ray scatter. Shielding area is ~ 30" (76cm) x 23.5" (60cm). Does not include ceiling suspension arm.</p>	A610-0215
	<p><b>Compass Contour and Frame Drape</b> Custom shaped drape for the Compass Contour and Frame Radiation Shield. Slips over the entire shield and is held in place with an elastic band. Sterile. Case of 20.</p>	C610-0256
	<p><b>IV Pole C000-1324</b> Stainless steel with four bottle hooks Range: 1070 -1660 mm Tube diameter: 25 / 18 mm <b>IV Pole Accessory Rail Clamp C000-1325 included</b></p>	A100-3285

## 100E Table Optional and Replacement Accessories:

Photo	Description <span style="color: red;">Red text = Metric accessory rails</span>	Catalog Number
	<p><b>Simple clamp</b> Accommodates a round accessory post. All metal design clamp mounts and locks anywhere on standard OR table rail. Autoclavable.</p>	C100-3046
	<p><b>Snap-on Blade Clamp</b> Accommodates a blade-style accessory post. Quickly snaps onto rail, push button to release clamp. Mounts and locks anywhere on standard OR table rail. Prevents false rail engagement. Autoclavable.</p>	C100-3043
	<p><b>Simple Blade Clamp</b> Accommodates blade-style accessory post. Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost-effective solution.</p>	C000-1115
	<p><b>Snap-on Socket</b> Accommodates a round accessory post. Quickly snaps onto rail and pull bottom lever to release. Mounts and locks anywhere on standard OR table rail. Lithotomy wheel for accessory angle adjustments.</p>	C000-0746
	<p><b>Restraint Straps</b> Set of two patient restraint straps with buckle locking mechanism and one pad for patient safety and comfort</p>	A100-3731

## 100E Shipping Crate Specifications

- Weight (approx.): 495 Lb (225kg)
- Length: 96" (2438mm)
- Width: 42" (1067mm)
- Height: 47" (1194mm)

## Pain Management / Spinal Imaging Table Specifications & Accessories

•Standard feature ◊Optional      NA Not available	<b>100E</b>	<b>100T</b>	<b>100RT</b>	<b>100RTL</b>
<b>Elevate (Height with 2" pad)</b>	27- 38.8" 68 - 101 cm	30 - 41.8" 76 - 106 cm	32.4 - 44.2" 82 - 112 cm	33.2 - 43.2" 84.3 - 109.7 cm
<b>Trendelenburg Tilt</b>	NA	±15°	±15°	±15°
<b>Lateral Tilt</b>	NA	NA	20°	20°
<b>Longitudinal Tabletop Motion</b>	NA	NA	NA	10"
<b>Standard Tabletop Size</b>	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm
<b>Optional Tabletop Size</b>	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out	22" x 80" 56 x 203 cm with facial cut-out
<b>Optional GI Tabletop Size</b>	28" x 80" 71 x 203 cm	28" x 80" 71 x 203 cm	28" x 80" 71 x 203 cm	28" X 80" 71 x 203 cm
<b>Imaging Area</b>	66" 167.6 cm	62" 157.5 cm	62" 157.5 cm	59" 149.9 cm
<b>Patient Capacity</b>	450 Lb 204 kg	450 Lb 204 kg	450 Lb 204 kg	450 Lb. 204 kg
<b>Table Weight</b>	290 Lb 131 kg	325 Lb 147 kg	380 Lb 172 kg	355 Lb. 161 kg
<b>Tabletop Pad (2" thick)</b>	•	•	•	•
<b>Restraint Straps</b>	•	•	•	•
<b>Facial Support Cushion</b>	•	•	•	•
<b>Pendant Hand Control</b>	•	•	•	•
<b>Foot Control</b>	◊	◊	◊	◊
<b>Armboard, rail mounted</b>	◊	◊	◊	◊
<b>Armboard, carbon fiber</b>	◊	◊	◊	◊
<b>VAB Armboard</b>	◊	◊	◊	◊
<b>Clamp-on Accessory Rail</b>	◊	◊	◊	◊
<b>Scatter-Guard Radiation Shield</b>	◊	◊	◊	◊
<b>Curve Spine Frame</b>	◊	◊	◊	◊



Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

**Image Diagnostics, Inc.**

310 Authority Drive  
Fitchburg, Massachusetts  
01420

For customer service:  
[sales@imagediagnostics.com](mailto:sales@imagediagnostics.com)  
Telephone: 978.829.0009



Image Diagnostics, Inc. 310 Authority Drive Fitchburg, MA, 01420 USA  
A veteran-owned company since 1986. Designed and manufactured in the U.S.  
Tel: 978-829-0009 Image Diagnostics, Inc. reserves the right to update product specifications without notification. Warranty information for IDI products is available upon request from Image Diagnostics, Inc.

