

Interconnect cable guide



WARNING!

Read the operator manual prior to use



Scan Here for
Operator Manual

idi
imagediagnostics.com

Included adaptors

J000-0292

DVI to HDMI Cable – 8in.



J000-0307

Display Port to HDMI Cable – 6in.



Order Replacements

1. Contact IDI at **978-829-0009**. Press 1 for Technical Support.

2. Have the **Serial Number** and **Part number** ready.

Business hours are 8 am—5 pm EST.

After-hours or weekend requests can be sent via email to techsupport@imagediagnostics.com. Please include the device's serial number, the required parts, the facility, name, phone number, and email address for prompt response when the business reopens the following business day.

Included Cables

J000-0346

Display Port to Display Port – 25ft.



J000-0293

HDMI Cable – 25ft.



J000-0236

Coaxial BNC Cable – 25ft.



J000-0186 (ilex only)

SVGA Cable – 20ft.



Connector Guide



WARNING!

Read the operator manual prior to use

Connector Types Provided

DVI-D




Digital

Primary Application

Most C-ARMS
Endoscopes
PC
Vitals
Ultrasound

VGA



Analog

Primary Application

PC
Vitals
Ultrasound

BNC



Digital and Analog

Primary Application

9900 - C-ARM
Endoscopes

HDMI



Digital

Primary Application

C-ARM w/adaptors
PC
Vitals

D-PORT



Digital

Primary Application

GE Elite C-ARM
ProVation/PACS PC