

Urology and GYN/Fertility Tables

100-UCPLUS 100US

Tables & Accessories



100UCPLUS Urology Table



100US GYN / Fertility Table



INTELLIGENT DESIGN ► TANGIBLE OUTCOMES

idi
imagediagnostics.com

100-UCPlus Urology/Multi-purpose Table



Shown with leg holders (#A100-2247).



Shown with tabletop extension (#A100-1400)

100UCPlus description:

Imaging table for urology and multi-purpose use. Urological accessories are easily removable so the table can be quickly configured for multi-purpose C-arm imaging or general patient positioning for non-urological procedures.

- 120VAC is standard
- 230VAC is optional and must be specified at time of order

100UCPlus #A100-1517 (120VAC with U.S. Accessory Rails)
A100-2639 (230VAC with metric accessory rails)

Includes the following:

- A100-1876: Stirrups Mounting Rail extension Kit
- X100-1818: tabletop pad, 28" (711 mm) x 46 ½" (1181 mm) x 2" (51mm) thick
- A100-1400: Tabletop extension, 24" (610 mm)W x 34" (864 mm) (Metric tables A100-2675)
- X100- 1817: Tabletop extension pad, 24" (610 mm) x 34" (864 mm) x 2" (51 mm) thick
- K000-0158: Pendant hand control
- A100-1540: Foot switch control
- U100-3025: Removable x-ray shield flap, for perineal end of table
- A100-3883: Patient restraint strap with pad
- A310-056: Arm board, rail mounted, with pad and strap (A310-062 for metric rails)
- C000-0595: Urology drain bag support hoop
- C000-0593: Disposable drain bag & hose assembly, sterile, box of 20.
- C000-0492: Disposable covers for foot switch, box of 50
- C000-0645: Disposable tabletop pad protector (qty of 1)
- C000-0612: 5 gallon collection bag (qty of 1)

100UCPlus Package Special #A100-2521 (120VAC)
A100-3180 (230VAC)




Includes the following:

- A100-1517: Aspect 100UCPlus Urology Table with standard accessories
- A100-2246: Two ComfortFin® 350 stirrups with clamps for US rails
- A310-056: One rail-mounted arm board
- C000-0593: Disposable drain bag & hose assembly, sterile, box of 20
- C100-3044: One ComfortFin stirrup storage cart

100UCPlus Package Special #A100-3107 (120VAC)
A100-3181 (230VAC)

Includes the following:

- A100-1517: Aspect 100UCPlus Urology Table with standard accessories
- A100-2247: Two ComfortFin 500 stirrups with clamps for US rails
- A310-056: One rail-mounted arm board
- C000-0593: Disposable drain bag & hose assembly, sterile, box of 20
- C100-3044: One ComfortFin stirrup storage cart

100UCPlus Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	IDI ComfortFin 350 Lithotomy Leg holders (Stirrups pair only) ComfortFin Lithotomy Leg holders are designed to improve infection control with durable, easy to clean silicone straps. Includes integrated lateral fin to help protect the head of the fibula and the peroneal nerve. 4" more working space with an increased range of motion. Patient weight capacity: 350 lbs (159 kg).	C100-3042
	Replacement Boot with Straps for IDI C100-3042 ComfortFin 350 Lithotomy Leg Holders	C100-3080
	Kit: IDI ComfortFin350 Lithotomy Leg holders C100-3042 with C100-3046 economy side-rail clamps	A100-2246 A100-3074
	IDI ComfortFin 500 Lithotomy Leg holders (Stirrups pair only) ComfortFin Lithotomy Leg holders are designed to improve infection control with durable, easy to clean silicone straps. Includes integrated lateral fin to help protect the head of the fibula and the peroneal nerve. Boots are sized for today's heavier patients. 4" (102mm) more working space with an increased range of motion. Patient weight capacity: 500 lbs. (227 kg).	C100-3045
	Replacement Boot with Straps for IDI C100-3045 ComfortFin 500 Lithotomy Leg Holders	C100-3079
	Kit: IDI ComfortFin 500 Lithotomy Leg holders C100-3045 with C100-3046 economy side-rail clamps.	A100-2247 A100-3075
	Lithotomy Leg Holder Cart Easy to push/pull with ergonomic handles. Stirrups are isolated to prevent damage. Clamp storage compartment. Designed to safely store ComfortFin stirrups when not in use.	C100-3044

100UCPlus Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Replacement Tabletop Pad 2" (51mm) thick foam pad. 28" (711mm) wide x 46 1/2" (1181mm) long. Replacement Extension Tabletop Pad 2" (51mm) thick foam pad. 28" (71mm) wide x 34" (864mm) long.	
	Physician Elbow Supports (pair) Removable	C000-0594
	Single use soft armband straps, package of 12	C000-1555
	Split Leg Holders For use with Aspect 100UC Plus Tables. Includes pads. Prone percutaneous, lateral or modified lateral nephrostomy. Includes two split leg holder clamps # C000-1115 or C000-1217 .	C000-1146 C000-1216
	Right replacement pad for Split Leg Holder C000-1146 or C000-1216 .	C000-1357
	Left replacement pad for Split Leg Holder C000-1146 or C000-1216 .	C000-1358
	Replacement Split Leg Holder Clamp For use with split leg holders. # C000-1146 and C000-1216	C000-1115 C000-1217
	Simple clamp Accommodates a round accessory post. All metal design clamp mounts and locks anywhere on standard OR table rail. Autoclavable.	C100-3046
	Snap-on Blade Clamp Accommodates a blade-style accessory post. Quickly snaps onto rail, push button to release clamp. Mounts and locks anywhere on standard OR table rail. Prevents false rail engagement. Autoclavable.	C100-3043

100UCPlus Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Simple Blade Clamp Accommodates blade-style accessory post. Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost-effective solution.	C000-1115
	Snap-on Socket Accommodates a round accessory post. Quickly snaps onto rail and pull bottom lever to release. Mounts and locks anywhere on standard OR table rail. Lithotomy wheel for accessory angle adjustments.	C000-0746
	Armboard Quick-release, rail mounted (requires Clamp-on Accessory Rail A100-1007 or A100-2379 for use on tabletop extension.	A310-056 A310-062
	Replacement pad for A310-056 or A310-062 Armboard	C100-3329
	Clamp-on Accessory Rail (For use on tabletop extension) Rail is 8 1/2" long x 3/8" thick x 1-1/8" high. Clamp-on Accessory Rail Metric rail is 10 mm x 25 mm x 216 mm	A100-2218 A100-2838
	Stirrups Mounting Rail extension Kit Included with table. For use with leg-holders (Leg-holders are sold separately).	A100-1876 A100-2676
	Hand Control for 100UCPlus and 100UC. Controls height up/down, trendelenburg/reverse, slide, cross table movements. (Included with table).	K000-0158
	Foot Switch Control for 100UCPlus. Included with table.	A100-1540

100UCPlus Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Compass Contour Radiation Shield Contoured ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhance eye, head and neck protection from x-ray scatter. Shielding area is ~30" (760 mm) x 23.5" (600 mm) with a 10" (254 mm) radius along the bottom edge of the leaded acrylic shield. Does not include ceiling suspension arm.	A610-0200
	Compass Frame Radiation Shield Rectangular ceiling suspended acrylic barrier radiation protection with a lead equivalence of 0.5mm. Enhance eye, head and neck protection from x-ray scatter. Shielding area is ~30" (760 mm) x 23.5" (600 mm). Does not include ceiling suspension arm.	A610-0215
	Compass Contour and Frame Drape Custom shaped drape for the Compass Contour and Frame Radiation Shield. Slips over the entire shield and is held in place with an elastic band. Sterile. Case of 20.	C610-0256
	IV Pole Stainless steel with four bottle hooks Range: 42" (1070 mm) / 65" (1660 mm) Tube diameter: .98" (25 mm) / 0.079" (18 mm) IV Pole Clamp C000-1325 is included.	A100-3285
	Vinyl Urology Table Covers Box of 20.	C000-0645
	Uro Drain Bags with Hose Disposable drain bag and hose assembly. Box of 20.	Sterile: C000-0593 Non-sterile: C000-1111
	Disposable Covers for Aspect 100UCPlus Foot Controls 24" (610 mm) x 20" (508 mm) Box of 50.	C000-0492

100UCPlus Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Uro Drain Support Hoop Collapsible & removable drain bag support hoop.	C000-0595
	Restraint Strap Patient restraint strap with hooks to attach to accessory rails on each side of the tabletop, a pad for patient safety and comfort and buckles to adjust to the patient girth	A100-3883
	Uro Collection Unit 5 gallon/18.9 liter capacity, Box of 10. Soft sided Vinyl container gently expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Built in zipper and handles allow for storage, closure of drain bag and hose assembly for enhanced infection control during transport and disposal.	C000-0612
	Uro Collection Unit 15 liter capacity, Box of 10. Flexible plastic sided container expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Integrated valve system allows user to shut off the flow of liquid from the drain bag in order to change collection units.	C000-1307
100 UCPlus Shipping Crate Specifications		
<ul style="list-style-type: none"> • Weight (approx.): 815 lbs. (370 kg) • Length: 96" (2438 mm) • Width: 42" (1067 mm) • Height: 47" (1194 mm) 		

Aspect 100-US Gyn/Fertility Clinic Table



Aspect 100US A100-1479 (120VAC with U.S. accessory rails)
A100-2674 (230VAC with metric accessory rails)

GYN/Fertility Mobile Imaging Table


An economical C-arm table for hysterosalpingograms

- 120VAC is standard.
- 230VAC is optional and must be specified at time of order

Includes the following:

- A100-758: Foot control
- X100-1687: Tabletop pad, 24" x 46" (610 mm x 1168 mm) 1" (25.4 mm) thick
- A100-3883: Patient restraint strap
- A100-2380: Leg support tabletop extension with pad
- A100-2248: Basic Knee Crutches [pair] with a pair of rotating side rail clamps [C000-0746]

100US Optional and Replacement Accessories:

Photo	Description	Red text = Metric accessory rails	Catalog Number
	Replacement Tabletop Pad 2" (51 mm) thick foam pad. 24" (610 mm) wide x 46" (1168 mm) long.		X100-1687
	Armboard Quick-release, rail mounted.		A310-056 A310-062
	Replacement pad for A310-056 or A310-062 Armboard		C100-3288
	Single use soft armband straps, package of 12		C000-1555
	Hand control for 100US		K000-0041

100US Optional and Replacement Accessories:		
Photo	Description Red text = Metric accessory rails	Catalog Number
	Basic Knee Crutches For patients up to 400 lbs. (180 kg) Basic pair with pads. Includes two C000-0746 rotating side rail clamps.	A100-2248 A100-2677
	Foot Control 100US (Included with table.)	A100-758
	IV Pole Stainless steel with four bottle hooks Range: 42" (1070 mm) / 65" (1660) mm Tube diameter: .98" (25 mm) / .079" (18 mm) IV Pole Clamp C000-1325 is included.	C000-1324
	Compass Contour Radiation Shield Contoured ceiling suspended acrylic barrier radiation protection with a lead equivalence of .5mm. Enhance eye, head and neck protection from x-ray scatter. Shielding area is ~30" (760 mm) x 23.5" (600 mm) with a 10" (254 mm) radius along the bottom edge of the leaded acrylic shield. Does not include ceiling suspension arm.	A610-0200
	Compass Frame Radiation Shield Rectangular ceiling suspended acrylic barrier radiation protection with a lead equivalence of 0.5mm. Enhance eye, head and neck protection from x-ray scatter. Shielding area is ~30" (760 mm) x 23.5" (600 mm). Does not include ceiling suspension arm.	A610-0215
	Compass Contour and Frame Drape Custom shaped drape for the Compass Contour and Frame Radiation Shield. Slips over the entire shield and is held in place with an elastic band. Sterile. Case of 20.	C610-0256

100US Optional and Replacement Accessories:			
Photo	Description	Red text = Metric accessory rails	Catalog Number
	Restraint Strap Patient restraint strap with hooks to attach to accessory rails on each side of the tabletop, a pad for patient safety and comfort and buckles to adjust to the patient girth		A100-3883
	Simple clamp Accommodates a round accessory post. All metal design clamp mounts and locks anywhere on standard OR table rail. Autoclavable.		C100-3046
	Snap-on Blade Clamp Accommodates a blade-style accessory post. Quickly snaps onto rail, push button to release clamp. Mounts and locks anywhere on standard OR table rail. Prevents false rail engagement. Autoclavable.		C100-3043
	Simple Blade Clamp Accommodates blade-style accessory post. Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost-effective solution.		C000-1115
	Snap-on Socket Accommodates a round accessory post. Quickly snaps onto rail and pull bottom lever to release. Mounts and locks anywhere on standard OR table rail. Lithotomy wheel for accessory angle adjustments.		C000-0746
	Uro Collection Unit 5 gallon/18.9 liter capacity capacity, Box of 10. Soft sided Vinyl container gently expands into a cube as it collects waste fluid when floor drain or suction is unavailable. Internal baffle system enables container to fill evenly as powder converts liquid to gel. Built in zipper and handles allow for storage, closure of drain bag and hose assembly for enhanced infection control during transport and disposal.		C000-0612

100US Optional and Replacement Accessories:			
Photo	Description	Red text = Metric accessory rails	Catalog Number
100US Shipping Crate Specifications			
<ul style="list-style-type: none"> Weight (approx.): 815 lbs. (370 kg) Length: 96" (2438 mm) Width: 42" (1067 mm) Height: 47" (1194 mm) 			

Aspect Urology & GYN Table Specifications & Accessories

	100 UC-Plus	100 US
Elevate (Height with pad)	33.3" – 43.1" (845 mm – 1094 mm)	27" – 39" (686 mm – 991 mm)
Trendelenburg Tilt	± 15°	NA
Transverse travel	± 3.5", 7" total (± 89 mm, 178 mm total)	NA
Longitudinal Tabletop Motion	10" (254 mm)	NA
Tabletop Size	28"W x 46"L (710 x 1170 mm)	24"W x 46"L (610 mm x 1168 mm)
Tabletop Extension Size	24"W x 34"L (610 x 860 mm)	NA
Imaging Area	Main tabletop: 25"W x 20"L (635 mm x 508 mm)	22"W x 32"L (559 mm x 813 mm)
Imaging Area with table extension	54" (1370 mm)	NA
Patient Weight Capacity	Upper body on main table: 525lb. (238 kg) Upper body on extension: 400 lbs (181kg) Load capacity on extension only: 200 lbs (90.7 kg)	450 lbs (204 kg)
Table Weight	490 lb. (222 kg)	375 lbs ((170 kg)
Tabletop Pad 2" (51mm) thick	•	•
Tabletop extension with pad	•	For leg support only
Hand control pendant	•	•
Footswitch control	•	◇
Patient restraint straps	•	•
Removeable drain bag support hoop	•	NA
Disposable drain bag & hose assembly, box of 20	•	NA
Removable radiation shield flap for perineal end of table	•	NA
Armboard, rail-mounted with pad and strap	•	◇
Split Leg Stirrups	◇	NA
Disposable covers for footswitch, box of 50	•	◇



Since 1986, IDI has been a manufacturer of innovative products for surgical, interventional and diagnostic imaging.

Our product focus is in electro-mechanical solutions for fluoroscopic and digital radiographic applications.

IDI is committed to superior product design, value and customer support in everything we do.

Image Diagnostics, Inc.
310 Authority Drive
Fitchburg, Massachusetts
01420

For customer service:
sales@imagediagnostics.com
Telephone: 978.829.0009

INTELLIGENT DESIGN ► TANGIBLE OUTCOMES



ORDER: sales@imagediagnostics.com
SUPPORT: techsupport@imagediagnostics.com
WEBSITE: <https://imagediagnostics.com>

+1-978-829-0009
Monday through Friday



This page intentionally blank

This page intentionally blank

This page intentionally blank

