



ComfortFin[®] Stirrups

Lithotomy positioning just got easier *and* safer!



ComfortFin 350 Stirrups
on IDI 100UC+ Urology Table



**ComfortFin
Stirrup Cart**



◀ Safer Prep

Medial to lateral foot cover enclosure keeps dangling straps away from the prep area

Safer Cleaning ▶

Silicone straps should be part of your comprehensive Bioburden Reduction Plan.

Silicon closure straps are much easier to clean and more durable than VELCRO® style straps.



◀ Better Pressure Protection

ComfortFin Stirrups have superior padding compared to standard stirrups.

The padding wraps the boot shell, minimizing fascia injury from hard plastic contact.

Better Sizing ▶

Unlike traditional stirrups, the ComfortFin 350 & 500 have boots sized to accommodate larger patients. Correct sizing reduces risk of pressure injuries.



Complete Urology Package



with ComfortFin™ Stirrups

PACKAGE INCLUDES:

- IDI's 100UC+ urology table and standard accessories*
- Two ComfortFin stirrups: (choose size)
 - ComfortFin 350 Stirrups (patients up to 350lbs)
 - OR
 - ComfortFin 500 Stirrups (patients up to 500lbs)
- Two basic side rail clamps (US rails)
- One rail-mounted quick-release armboard
- One box of twenty Uro Drain bags (sterile) with hose/disposable drain bag assembly
- One ComfortFin stirrup cart
 - Designed to securely store ComfortFin Stirrups

IDI's 100UC+

The industry's best-selling
Urology table!



ComfortFin 350 or
500 Stirrups

A100-2246- ComfortFin 350

A100-2247- ComfortFin 500



ComfortFin Storage
Cart

C100-3044



Quick-release
armboard

A310-056



Uro Drain Bags
(box of 20)

C000-0593



Quick-release
Blade Clamp

C100-3046

ORDERING INFO:

ComfortFin 350 package

Part #: A100-4114

ComfortFin 500 package

Part #: A100-4115

SCAN

visit IDI's
ComfortFin
webpage



Standard Accessories included with 100UC+

- 2 inch thick tabletop pad
- Aluminum boarding extension
- Hand control pendant and wireless footswitch control
- Removable drain bag support hoop
- Disposable drain bag and hose, box of 20
- Patient restraint straps, qty 2
- Removable radiation shield flap for perineal end of table
- Armboard, rail-mounted with pad and strap
- Disposable footswitch covers, box of 20



Image Diagnostics, Inc. 310 Authority Drive Fitchburg, MA, 01420 USA
A veteran-owned company since 1986. Tel: 978-829-0009 Fax: 978-829-0027
imagediagnostics.com ComfortFin® is a registered trademark of Image Diagnostics Inc.
Velcro® is a registered trademark of Velcro BVBA ¹Tested 27,000 cycles with no failure.
Hook and loop loses 50% of holding strength at 2,500 cycles. L100-3116C



imagediagnostics.com